















**CABOS ECG E SPO2 PARA EQUIPAMENTO HOSPITALAR**


|  |   |   |  |
|--|---|---|--|
| <p><b>ML4491</b><br/>Lohmeier Multi-parameter Trunk Cable, (IBP*2, Temperature*2, 5-lead ECG , Spo2 ) 98.060.120/A<br/>Lohmeier Multi-parameter Trunk Cable 98.060.120/A<br/>Lohmeier Multi-parameter Trunk Visicon M010, Cable 98060120/A</p>   |    | <p>Lohmeier Multi-parameter ECG Trunk Cable 5Ld Multi-Cable, 9ft IEC/AAMI, 27pin plug&gt;Euro 5Ld IBP*2, Temperature*2, ECG, Spo2 (discharge inside) Compatible Models No. Lohmeier D-81241, P/N.: Lohmeier 98.060.120/A, (98060120) Monitor B. Braun M010 Ref.98060120</p> | <p>Use <b>ML4919</b><br/>LOHMEIER ECG leadwire set , 5-Lead , Grabber , IEC, 40IN &amp; 2pin</p>   |
| <p><b>ML4919</b><br/>LOHMEIER ECG leadwire set , 5-Lead , Grabber , IEC, 40IN &amp; 2pin To fit Lohmeier Multi-link trunk cable ML4491 ECG parameter</p>   |    |    | <p>To fit Lohmeier Multi-link trunk cable <b>ML4491</b> ECG parameter</p>  |
| <p><b>ML4492</b><br/><b>5590539 DRAGER</b><br/>Multi-Parameter NEOMED Trunk Cable (temp*2, ECG, SPO2, FIO2 ) - Siemens 5590539<br/>Siemens drager 5590539</p>  |    | <p>DRAGER 5590539 Multi-Parameter NEOMED Trunk Cable Multiparameter cable to monitor neonatal 3-lead ECG, impedance respiration, SpO2, two temperatures (Infinity modular monitors only), and FiO2 (Infinity modular monitors only).</p>                                    | <p>Not for use during Electrosurgery Including: - Pod mount.</p>   |
| <p><b>ML4493</b><br/><b>3368391</b>, Cabo ECG c/5 derivações p/monitor EAGLE Mod. 3100 e (P/monitor Siemens SC 7000, SC 6000, SC 7000, SC 9000)( p/ECG Siemens Sirecust 1260). Siemens SC9000 XL Multi-Link Trunk cable (Spo2, Temperature, Ecg). Use with Siemens Instrument: SC6000, SC6002XL, SC7000, SC9000, Siemens Art. 3368391, 3368391 B, 3368391 0B. 3368391 SC6000, SC6002XL, SC7000, SC9000, Siemens 3368391 Siemens drager 3368391</p> |   |   | <p>Siemens SC9000 XL Multi-Link Trunk cable (Spo2, Temperature, ECG (SC 7000) ICU/CCU. Marque : DRAGER Multimed 5 1,5m ref. 59 50 169 E5300U, 5950169E5300 U, for monitor Siemens SC 7000)</p> |
| <p><b>ML4519</b><br/><b>5956433</b>, Siemens ECG leadwire set, 3-Lead Snap, IEC, 40IN &amp; 2pin. Cabo p/monitor Siemens drager ref.5956433. Siemens drager 5956433. Siemens 5956433</p>   |  | <p>Siemens Instrument SC6000, SC6002XL, SC7000, SC9000,</p>    | <p>Cabo ECG Précordial (1M) de 3 cabos clip 5956433</p>  |
| <p><b>ML1588</b><br/><b>7486942</b>, Siemens ECG Leadwire set, 3-Lead Snap, IEC, 40IN &amp; 2pin. Cabo p/monitor Siemens drager ref.5956433. Cabo p/ Monitor Siemens / Drager ref. 7486942 Siemens, 3 derivações. Cabo p/monitor Siemens drager ref.59564333 (AK3) Siemens 7486942</p>   |  | <p><b>ML1588</b><br/>Siemens Series Monitors 6000/6002/6802/7000/9000 and Dräger monitors gamma / GammaXL / Delta (Multimed pod system for Series SC 6000, SC 7000, SC 9000). Siemens drager 7486942</p>  |  |
| <p><b>ML3302</b><br/>5956466E530U Siemens ECG Leadwire set 5-Lead, Grabber, IEC, 40IN &amp; 2pin, 5956466E530U. Cabo Précordial Monitor PM-7000</p>   |  | <p><b>ML3302</b><br/>Siemens Instrument :SC6000, SC6002XL, SC7000, SC9000,</p>   | <p>Cabo ECG Précordial (1M) de 5 cabos clip: 5956466. Siemens 5956466</p>  |




|  |   |  |   |
|--|---|--|---|
| <p><b>ML3303</b><br/>Datex-Ohmeda S/5, 10-lead Trunk cable , Green connector , 12pin female connector<br/>Datex-Ohmeda GE<br/>Use with leadwire <b>ML3301</b> and <b>ML3301C</b></p>   |    |   | <p><b>ML3301</b><br/>ECG leadwire set , <b>5-Lead , Snap, IEC</b>, Connect trunk cable <b>RA, LA, RL, LL, V</b><br/><b>ML3301C</b><br/>ECG leadwire set, 5-Lead Snap, IEC, 40IN Connect trunk <b>Cable C2-C6.</b></p> |
| <p><b>ML4474</b><br/>Cable GE Multilink M1099319, 12-ld ECG Trunk Cable, <b>IEC</b>, 3m/10ft (3M)</p>  <p><b>C2-C6</b>      <b>N-C1</b></p>  |    |  <p>See Also <b>5-Lead Snap MI1146</b> and <b>ML1147</b></p>                                      | <p><b>ML1146G</b><br/>Leadwire set, 5-Lead, <b>N-C1 grabber, IEC</b>, 29in (76cm)</p> <p><b>ML1147G</b><br/>Leadwire set, 5-Lead, <b>C2-C6, grabber, IEC</b>, 29in (76cm)</p>   |
| <p><b>ML5437</b><br/>GE Eagle, Solar, Dash ,Tram , MAC-lab cath lab, Datex Ohmeda S/5 FM, GE One-piece EKG CABLE and leadwires, <b>IEC, Banana</b>, 11pin male</p>   |     | <p>Cabo GE Eagle, Solar, Dash ,Tram , MAC-lab cath lab, S/5 FM, GE One-piece EKG CABLE and leadwires , <b>IEC, Banana</b> , 11pin male</p>   |   |
| <p><b>ML5431</b><br/>416035-002, GE Eagle, Solar, Dash, Tram, MAC-lab cath lab, Datex Ohmeda S/5 FM, Multi-Link 12-Lead ECG Standard <b>Cable, IEC</b>, 12 ft.</p> <p><b>ML5432</b><br/>416035-001, GE Eagle, Solar, Dash ,Tram , MAC-lab cath lab, Datex Ohmeda S/5 FM, Multi-Link 12-Lead ECG Standard <b>Cable, AHA</b>, 12 ft.</p> |    <p><b>C2-C6</b>      <b>N-C1</b></p> |  <p>See Also <b>5-Lead Grabber MI1146G</b> and <b>ML1147G</b></p>                                | <p><b>ML1146</b><br/>Leadwire set, 5-Lead, <b>N-C1 Snap, IEC</b>, 29in (76cm)</p> <p><b>ML1147</b><br/>Leadwire set, 5-Lead, <b>C2-C6, Snap, IEC</b>, 29in (76cm)</p>   |
| <p><b>ML2111</b><br/>M1500A, M1510A, M3046A, HP 3-Lead monitoring IEC trunk cable, 9ft &amp; 12pin.</p>   |     |  <p><b>ML2112</b><br/><b>ML2112A</b><br/><b>ML2112AF</b></p>                                    | <p><b>ML2112</b><br/>Cable M1613A, HP ECG safty Lead set 3-Lead, IEC, 3ft, AA 2PIN/Grabber.</p>   |
| <p><b>ML6707</b><br/>M1500A, M1510A, M3046A, HP 3-Lead monitoring <b>AHA</b> cable, 9ft &amp; 12pin.<br/><b>One piece 3-lead cable, grabber</b></p>  |     | <p>Philips Medical System: 78352C, 78354C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A. Philips 989803143171. One piece 3-lead cable, grabber</p> |   |
| <p><b>ML2002</b><br/>M1520A, M1530A, HP 5-Lead monitoring IEC trunk cable 9ft &amp; 12pin. M1510A, M3046A.</p>    |     |   | <p><b>ML1002</b><br/>Précordial cable for M1520A, M1530A M1633A</p>   |
| <p><b>RL1115</b><br/>Philips MP70, MP70, M1668A Cable 5 Lead ECG Trunk, AAM/IEC, label included. Round 12PIN 2.7m. HP M1668A 5-lead Trunk cable. Usa Cabo ML1008.</p>  |                                        |   | <p><b>ML1008</b><br/>( M1971A ) HP M1971A 5-lead leadwires, IEC, Grabber, Para Monitor Philips MP70. Usa o</p>  |




|  |   |  |  |
|--|---|--|--|
| <p><b>RL1117</b><br/>Philips MP70, HP M1669A, HP-2385, HP 2385, 3-lead Trunk cable, AAMI &amp; IEC label included. Round 12PIN</p>   |                                        |   | <p><b>RL1119</b><br/>HP M1974A 3-lead leadwires, IEC, Snap.</p> <p><b>RL4499</b><br/>HP M1673A 3-lead leadwires, AHA, Snap</p>             |
| <p><b>HP-2385</b><br/>3-lead ECG patient trunk cable with 12-pin connector, AAMI and IEC, OEM P/N: M1669A, ECG Patient Cable Sets (3 Leads) can be just used with Philips monitoring like 78352C, 78354C, 78834C, M1001A/B, M1002A/B, M3001A, M3002A, M8105A, M3535A, M3536A, OEM P/N: M1949A + M1968A + M1976A.</p>   |                                        |   | <p>HP3-90P. (ICU grabber, AAMI, OEM P/N: M1671A) Multi-Measurement Modules M3000A, M3001A</p>  |
| <p><b>HP-21085</b><br/>10-Lead ECG Patient Cable Sets (5 + 5 Leads) HP-21085, HPB5-90P, HPA5-90P. Unimed 10-Lead ECG Patient Cable Sets (5 + 5 Leads) can be just used with Philips monitoring like 78352C, 78354C, 78834C, M1001A/B, M1002A/B, M3001A, M3002A, M8105A, M3535A/M3536A, OEM P/N: M1949A + M1968A + M1976A.<br/><b>Veja n/ref. ML112, ML9114</b><br/><b>See Our ref. ML112, ML9114</b></p> |                                        | <p>Unimed 10-Lead ECG Patient Cable Sets (5 + 5 Leads)</p>  <p><b>Veja n/ref. ML9113, ML9114</b><br/><b>See Our ref. ML112, ML9123</b></p>                        | <p>Unimed 10-Lead ECG Patient Cable Sets (5 + 5 Leads), are designed to be used with Philips New ECG Monitoring Cable and Lead System.</p> |
| <p><b>HPB5-90P.</b><br/>(ICU grabber, limb, AAMI, OEM P/N: M1968A. 78352C, 78354C, 78834C Legacy HP Patient Monitors<br/><b>Veja n/ref. ML9113, ML9114</b><br/><b>See Our ref. ML112, ML9123</b></p>   |                                      | <p>Unimed 10-Lead ECG Patient Cable Sets (5 + 5 Leads) can be just used with Philips monitoring like 78352C, 78354C, 78834C, M1001A/B, M1002A/B, M3001A, M3002A, M8105A, M3535A/M3536A, OEM P/N: M1949A + M1968A + M1976A.</p>                       |  |
| <p><b>ML1101</b><br/><b>DISPOSABLE HAND SWITCH PENCILS</b><br/>product name: Surgical apparatus<br/>clicks: 87 times<br/>Length of handle(cm):160mm<br/>Length of blade(cm):40mm<br/>Diameter of blade:2.3mm<br/>Length of wire(m):2.8m<br/>Sterilized by EO</p>   |  <p><b>ML1101</b> <b>ML1122</b></p> | <p><b>ML1122</b><br/><b>REUSABLE HAND SWITCH PENCILS</b><br/>Cable length: 320cm<br/>Ø shaft: 2,38mm<br/>Connector Valleylab<br/><b>Sterilization: 100</b></p>   |   |
| <p><b>ML2209</b><br/>Delfi Monitor 72303 ECG Cable. Delfi-Monitor 72303 ECG Cable (with 4-lead, IEC, CLIP)</p>   |                                      | <p>ECG cable (with 4-lead, IEC, CLIP) 72303 from Defi-Monitor ECO1. Primedic Defi-Monitor DM10/DM30 <b>4-lead ECG EKG Cable with leadwires for Medical equipment.</b> Primedic DM10, DM30 <b>4-lead ECG Cable</b> with leadwires, clip, ecg lead</p> | <p>Cabo ecg com 4 electrodos (leads) para desfibrilhador Primedic, Defi-Monitor ECO 1</p>  |
| <p><b>LM4422</b><br/>Cable Kenz Delta 3 plus Kanz EKG cable, IEC, DDK-16PIN &gt;Banana 4.0</p>   |                                      | <p>For KANZ 103. KANZ 104 KANZ 106. Cardioline Delta 1 plus, Cardioline Delta 3 plus, Cardioline Delta 60 plus. PC-104</p>   | <p>KANZ 103.104 106. Cardioline Delta 1 plus, Delta 3 plus, Delta 60 plus</p>  |




















|   |   |   |  |
|---|---|---|--|
| <p><b>ML6717</b><br/>Cabo ECG FUKUDA CP-100A<br/>para Eletrocardiógrafo Fukuda<br/>7110<br/>Fukuda 7110 ECG cable</p>   |     |  <p><b>ML2011</b></p>  | <p><b>ML2011</b><br/>Selector<br/>derivações</p>  |
| <p>Conector para ECG CP-100ª<br/>ECG 7110, Ref. <b>RC884422</b></p>    | <p>Conector para Cabo ECG CP-100ª, ECG 7110,<br/>Ref. <b>RC880099</b></p>                 | <p>Conector para ECG CP-100ª<br/>ECG 7110, Ref. <b>RC884422</b></p>                    | <p>For other<br/>material<br/>please<br/>contact us</p>  |
| <p><b>ML6055</b><br/>Cabo ECG FUKUDA CP-100J<br/>Eletrocardiógrafo Fukuda FK-12<br/>Fukuda FK-12 Cable ECG</p>  |    | <p>Cable for ECG Fukuda CP-<br/>100J/FK-12</p>  |  |
| <p><b>ML6591</b><br/>Cabo ECG Fukuda CP-50A<br/>MONOBLOC CABLES<br/>ECG cable with 5 leads banana<br/>Diameter connector 17mm<br/>Fukuda ME cable ,IEC, Round<br/>6pin &gt; 5 Leads Banana 4.0</p>          |     |   | <p>Cabo ECG<br/>Fukuda<br/>CP-50A<br/>6pin &gt;<br/>5-Lead<br/>Banana 4.0</p>  |
| <p><b>ML3173</b><br/>Fukuda Denshi ECG Cable,<br/>DSL-8001, DSL-8001R.<br/>Dynascope DS-8100,</p>   |     | <p>Dynascope DS-8200<br/>Dynascope DS-8500</p>  |  |
| <p><b>ML4422</b><br/>Galvanómetro para ECG<br/>Fukuda 7110. Parts<br/>Galvanometer for ECG 1 canal<br/>Fukuda Denshi FJC-7110<br/>Temos mais material de<br/>reparação</p>                                  |    |  <p><b>ML4429</b></p>  | <p><b>ML4429</b><br/>Botão Button<br/>derivações<br/>para ECG<br/>Fukuda<br/>7110. Parts</p>   |
| <p><b>ML4430</b><br/>FICHA BSS, Wandsworth Brand<br/>30A 250V B.S. 546 (BS546) (C653).<br/>30AMP BS 546 Roud Pin plug.<br/>Wandsworth 30A.250V B.S. 546<br/>Diameter 95x82,5mm</p>                          |    | <p><b>ML4430</b><br/>Wandsworth<br/>Plug 30A 250V B.S. 546<br/>(BS546) (C653). 30AMP<br/>BS 546 Roud Pin plug<br/>Diameter 95x82,5mm</p>                                  | <p>Wandsworth<br/>30A.250V B.S.<br/>546</p>  |
| <p><b>ML4812</b><br/>FICHA BSS, Wandsworth Brand<br/>Prod. 3034/W, Batch nº FZ5425,<br/>F25425, 15 AMP BSS<br/>Round Pin Plug, Ficha BSS 15A-FZ5425<br/>15 AMP BSS Wandsworth<br/>Wandsworth 15 AMP BSS</p> |    |    | <p><b>ML4812</b><br/>Wandsworth<br/>Plug 3034/W<br/>Batch nº<br/>FZ5425<br/>F25425.<br/>Round Pin Plug</p>                           |




|   |   |   |   |
|---|---|---|---|
| <p><b>ML1005</b><br/>Cabo ECG H.P. para Eletrocardiógrafo H.P. 1500 /1550A<br/>Cable for ECG H.P. 1500/1550A</p>  |    | <p>Stylus ECG H.P. 1500 /1550A</p>  <p><b>ML1003</b></p>   | <p><b>ML1003</b><br/>Estilete para ECG H.P. 1500 /1550A</p> |
| <p><b>ML9976</b><br/>Cabo para ecógrafo<br/>Ficha ligação equipamento 18 pins<br/>Ficha ligação sensor 16 pins<br/>Interface cable<br/>Device Connector 18 pin<br/>Sensor Connector 16 pin</p>  |   |   <p>Sensor Connector</p> | <p>Cabo para ecógrafia</p>                                  |
| <p><b>HT2793A</b><br/>D: Din lead wire ( Acc to yoke)<br/>5: 5-lead (3= 3 lead)<br/>90: 90cm(wire length)<br/>S: snap end(P=grabber end)<br/>I: IEC code (Others see below)</p>   |    |   |   |
| <p><b>ML5582</b><br/>Cables REF.70831 (60cm) ,<br/>70834 (1,20m) Grabber</p>  |   |   | <p>Cables<br/>REF.70852 (60cm)<br/>70855(1,20m)</p>         |
| <p><b>ML4555</b><br/>Cables REF.70811 (60cm)<br/>70814 (1,20m) Grabber</p>  |    |    | <p>Cables<br/>REF.70812 (60cm)<br/>7815(1,20m)</p>          |
| <p><b>ML1022</b><br/>Cabo précordial SNAP com 64cm comprimento, Diâmetro terminal 2mm</p>   |    |    |   |
| <p><b>ML4923</b><br/>ECG cable 3 leads for monitor DATEX OHMEDA</p>   |    |    | <p>Cabo ECG cable 3 leads for monitor DATEX OHMEDA</p>      |
| <p><b>ML9033</b><br/>Hellige 5-Ld ECG Trunk Cable. Spec:AHA/IEC, Round 10J&gt;ASP 5LD, 1K Resistor . Hellige, Servomed 4, Generation, SCP 851, SCP 852, SCP 844, SCS 804,, CardioServ, VICOM-sm</p>   |    |   |   |
| <p><b>ML1944</b><br/>S&amp;W (Artema) 5-Ld ECG trunk cable. Spec:AHA/IEC, Round 10J&gt;LL 5LD . S&amp;W (Artema) instrument: 8000 System, Diascope1, 2, Neoscope, DMS 700, defibrillator, Diascope-NT, Athena Systems, Cardio-Aid, DS Traveller, DSNT, DS1, DS2, MCX/LS</p> |    |   |   |
| <p><b>ML7811</b><br/>Cable 412682-003 GE Marquette Multi-link Leadwire set 3-Lead, Grabber, IEC2 9in (76cm) &amp; VS 2Pin. Cable 412682003 GE Marquette Multi-link Leadwire set 3-Lead</p>  |    |    |   |



|  |   |  |   |
|--|---|--|---|
| <p><b>ML7807</b><br/>Cable 545920, GE DATEX, OHMEDA, Integral process 71324. MONOBLOC, <b>3 LEADS CLIP</b>. ECG RESPIRATION, Old generation S5, CARDIOCAP - AS3 - CS3, 10 FEMALE PINS. ECG/Respiration cable, IEC color coding (545920),</p>           |   | <p>Specification: Three-lead ECG/respiration cable with clips, 3.5 m/11 ft., IEC color coding, green connector, Note: Use with Cardiocap and Cardiocap II.</p>   |   |
| <p><b>ML8511</b><br/>Cable 9030-01-49318-5. Diascope. Athena, SW 800, 72854, MONOBLOC <b>3 LEADS CLIP</b>. Artema SW SERIE 8000-ATHENA DIASCOPE 1-2 ATHENA DIASCOPE 2 S&amp;W 10 FEMALE CONTACTS. Cable 13179. <b>With 3 Leads Snap see ML8533</b></p> |    | <p>Artema SW, SERIE 8000-ATHENA, Diascope 1-2, 10 Female CONTACTS. Siemens Megacart PK200 (Siemens Elema AB) Cabo ECG p/Monitor Artema (Athena, SW) Transvalet 9030.01 de 3 terminais. (9030-01-49318-5)<br/><b>With 5 Leads Clip see ML6010</b><br/><b>With 5 Leads Snap see ML6033</b></p> | <p><b>72854</b><br/>3 Précordiais <b>(9030-01-49318-5)</b><br/>Diascope Athena, S&amp;W 800 ATHENA DIASCOPE 1 S&amp;W</p> |
| <p><b>ML5855</b><br/>EMBASE PEDIATRIQUE GAMME 8000-ATHENA Diascope 1-2, 10 CONTACTS FEMELLESML5852 ARTEMA SW, TRUNK CABLE PEDIATRIC. GAMME 8000-ATHENA DIASCOPE 1-2, 10 FEMALE PINS</p>  |    |   | <p><b>ML5612</b><br/>Terminais mola/ clips, Conjunto 5 Terminais 5 voies monitor Neotrak 502 70812</p>                    |
| <p><b>ML5852</b><br/>ARTEMA SW, 9030-10-49323-5, TRUNK CABLE PEDIATRIC. GAMME 8000-ATHENA DIASCOPE 1-2, 10 FEMALE PINS EMBASE PEDIATRIQUE GAMME 8000-ATHENA DIASCOPE 1-2, 10 FEMALE PINS</p>   |    |   | <p><b>ML4553</b><br/>Terminais mola/ clips, Conjunto 3 Terminais 3 voies monitor Neotrak 502 70812</p>                    |
| <p><b>ML2143</b><br/>72000, 545300, 545905 3 LEADS CLIP GE, DATEX OHMEDA. EMBASE, 3 VOIES SANS PROLONGATEURS. LIGHT-F-ML1 CARDIOCAP - AS3 - CS3 10 CONTACTS FEMELLES</p>   |  |   | <p><b>LM8822</b><br/>(1m) or <b>ML6655</b><br/>(1,5m)</p>   |
| <p><b>RH3304</b><br/>72002 (545301 or 545909) 5 LEADS MAIN CABLE WITHOUT LEADWIRE, DATEX, CARDIOCAP- AS 3 - CS 3</p>   |  |   | <p><b>ML2391</b><br/>Cabo ECG precordial de 5 cabos clip (Référence Constructeur 5956466)</p>                             |
| <p><b>ML2433</b><br/>Cable 71270 ECG (3 voies), for monitor Corometrics Neotrak 502. Cable ECG (3 voies) For monitor Corometrics NeoTrak 502</p>   |  |   | <p><b>LM8822</b><br/>(1m) or <b>ML6655</b><br/>(1,5m)<br/>(3 voies monitor Neotrak 502)</p>                               |
| <p><b>ML1906</b><br/>Mindray 040-000545-00 cable MR6701. MR6702 MR6905 Mindray BeneHeart Defibrillator Pads Cable With 50 ohm Test Load. Mindray cable 040-000545-00 cable</p>   |  |   | <p>Cabo de teste Mindray MR6701 com carga de 50 ohm 040-000545-00.</p>  |
| <p><b>ML5091</b><br/>Biolight M7000, M9500 <b>5-lead ECG</b> cable with leadwire, 12pin, IEC, Snap, one piece ECG cable with leadwire</p>  |  | <p><b>ML5092</b><br/>Biolight M7000, M9500 <b>3-lead ECG</b> cable with leadwire, 12pin, IEC, Snap, one piece ECG cable with leadwire</p>  |                                      |

















|   |   |  |   |
|---|---|--|---|
| <p>Cabo ECG p/monitor EAGLE Mod. 3100. AS73524MQ<br/>       Stammkabel kompatibel zu Marquette Hellige. 5-adrig für Solar, Dash, Eagle, Länge 2,75 m trunk cable , Marquette Hellige, 5 leads, for Dash, Eagle, Solar, Length: 2,75 m. AS414556-003<br/>       Patientenleitung . Original Marquette Hellige, 5 adrig mit Klammer, Länge: 74 cm.</p>  |    |   | <p>Marquette Hellige patientleads 5 leads with clamp, length: 107 cm.</p>   |
| <p><b>ML7981</b><br/>       (5-Lead) 2017003-003<br/>       GE-Marquette 5-LEAD ECG cable, IEC, 12ft, 11PIN, with 5-lead leadwire, Without resistance. Cabo ECG para Dash 2000 Marquett. de 3 derivações (reference for 3 leads is 72465 and reference for 5 leads is 72466). GE MARQUETTE HELLIGE EAGLE 3000 4000 MODULES TRAM DASH 2000 3000 4000. Type: DASH2000. I.P. Reference: 72465. MAIN CABLE 3 LEADS WITHOUT LEADWIRES. Old Manufacturer Reference: 407347002. EAGLE 3000/4000. 2017003-003. This trunk cable works with the following ECG monitors: Dash 1000, Dash Pro 2000, Dash Pro 3000, Dash Pro 4000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000, Solar 3000, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape. <b>ML7981</b>, Multilink E9002ZC Cabo ECG para Dash 5000</p> |  <p><b>ML7981</b></p>  <p><b>ML7077</b><br/>       Multi-link Leadwire set 5,-Lead 414556-003</p> | <p><b>ML7703</b><br/>       GE-Marquette 5-LEAD ECG cable ,IEC , 12ft,11PIN, with 3-lead leadwire</p>  <p><b>ML7811</b><br/>       Multi-link Leadwire set 3-Lead <b>ML7812: 412682-004</b></p> | <p>set 3-Lead</p>  <p>set 5-Lead</p>  |
| <p><b>ML7808</b><br/>       GE Medical, Dash Pro4000, DASH PRO 3000, DASH 2000,1000, Eagle system, MAC-Lab, Solar 9500, Solar 8000. One piece ECG cable with leadwire, 5 lead, IEC, Grabber. RC14426K, GE 2586P-1 Monitor Portátil Dash 2500 Cabo ECG para Dash 5000 (Direct Connect )</p>  |     |   | <p><b>ML9099</b><br/>       2021141-002 ECG cable with leadwires IEC, 3-Lead, Grabber Solar 3000, Tram 100.,200 , 300 , Tram 400, 500, 600. 2386P-I DASH 2500 Direct Connect</p>                            |
| <p><b>ML22011</b><br/>       JC-906P, Nihon Kohden JC-906P ECG trunk cable , 6-Lead , 12pin Plug to 8pin socket</p>   |     |   | <p>Nihon Kohden JC906P</p>  |
| <p><b>ML22044</b><br/>       BR-903P, Nihon Kohden BR-903P leadwires , 3-Lead, Grabber , IEC ,8pin socket</p>   |    | <p>BSM2300 Lifescope I/L BSM-4100 Lifescope P BSM-5100 Lifescope A BSM-9510 Lifescope M BSM-9800 Lifescope S BSM-1500 Lifescope N</p>  |    |
| <p><b>ML7774</b><br/>       M1735A, 3-Lead monitoring IEC one piece cable with leadwires, 10ft, 8pin.<br/>       Philips Medical, M1722A/B , M1723A/B, M2475B, M3500B, 78619A, 78660A/78660B, Defibrillator Cable to HP CODEMASTER 100 (M2475B)</p>   |     | <p><b>ML7770</b><br/>       M1736A, 5-Lead monitoring IEC one piece cable with leadwires, 10ft, 8pin.<br/>       Philips Medical, M1722A/B , M1723A/B, M2475B, M3500B, 78619A, 78660A/78660B, Defibrillator HP CODEMASTER 100</p>  |    |
| <p><b>ML4006</b><br/>       Nihon Kohden one piece ECG cable with leadwire6,-lead,IEC,Grabber, 12pin (Direct Connect )</p>  |     | <p>BSM2300 Lifescope I/L, BSM-4100 Lifescope P, BSM-5100 Lifescope A , BSM-9510 Lifescope M ,BSM-9800 Lifescope S , BSM-1500 Lifescope N</p>   |   |




|   |   |   |  |
|---|---|---|--|
| <p><b>ML4423</b><br/>Nihon Kohden one piece ECG cable with leadwire3,-lead,IEC,Grabber, 12pin<br/>(Direct Connect )</p>   |                                        | <p>BSM2300 Lifescope I/L<br/>BSM-4100 Lifescope P<br/>BSM-5100 Lifescope A<br/>BSM-9510 Lifescope M<br/>BSM-9800 Lifescope S<br/>BSM-1500 Lifescope N</p>   |  |
| <p><b>IP72809</b><br/>SIRECUST 600 - SIRECUST 732 (p/ECG Siemens Sirecust 1260) IP 72809, EMBASE 5 VOIES DRAEGER SIEMENS. SIRECUST 600 - SIRECUST 732 ( p/ECG Siemens Sirecust 1260). Sirecust 341-342 R , Series 200-400-600-700-900-1200, 10 Male Contacts</p>                                  |                                        |    |                                   |
| <p><b>IP72815</b><br/>SIRECUST 600 - SIRECUST 732 ( p/ECG Siemens Sirecust 1260) IP 72815, EMBASE 5 VOIES DRAEGER SIEMENS. SIRECUST 600 - SIRECUST 732 ( p/ECG Siemens Sirecust 1260). Sirecust 341-342 R , Series 200-400-600-700-900-1200, 10 Male Contacts</p>                                 |                                        |    | <p>Direct<br/>Connect</p>  |
| <p><b>ML8711</b><br/>Cables for ECG BURDICK ECLIPSE 8, Eclipse 850 (71125, ECLIPSE 4, ECLIPSE 8. BURDICK (QUINTON ) :Eclipse 4, Eclipse 400, Eclipse 8, Eclipse 800, Eclipse LE II,</p>   |                                        | <p><b>Burdick / QUINTON</b><br/>Eclipse 4, Eclipse 400, Eclipse 8, Eclipse 800, Eclipse LE II, Eclipse 850, Eclipse 850I, Eclipse Plus, Eclipse LE II, ATRIA, ATRIA PLUS, ATRIA 3000</p>                  |  |
| <p><b>ML5668</b><br/>Cabo ECG Drager infinity One piece ECG cable with 3-lead , AHA <b>Snap</b> leadwires .Round 7pin</p>   |  <p><b>ML5668</b> <b>ML5669</b></p>   | <p><b>ML5669</b><br/>Cabo ECG Drager infinity One piece ECG cable, 3-Lead , <b>Grabber</b>, AHA,12ft, Round 7pin</p>  |                                 |
| <p><b>ML5682</b><br/>Kontron 7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, Kolormon, Mcromon, Minimon, Supermon, Vitalmon. Vitalmon. One-piece ECG cable with leadwires , <b>3-Lead, Snap</b> IEC, 12pin</p>  |  <p><b>ML5682</b> <b>ML5683</b></p> | <p><b>ML5683</b><br/>Kontron 7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, Kolormon, Mcromon, Minimon, Supermon, Vitalmon. One-piece ECG cable with leadwires , <b>5-Lead Grabber</b>, IEC, 12pin</p> |                                 |
| <p><b>ML1948</b><br/>Biolight Q5 5 lead Ecg Cable, IEC, clip, grabber green 12pin Cabo ECG c/5 terminais p/monitor Biolight BLT Q5.</p>   |                                     |    |                                 |
| <p><b>ML4506</b><br/>Siemens / Hellige EKG cable, IEC, DB15M Banana 4.0. Cabo Siemens Megacart (Cabo Completo Monobloco 10 LEADS BANANA) Siemens Megacart (PK200) Elema AB, Hellige, Custo-med, Cambridge, DEGO ECG cable Cardio Control Perfect MD. <b>Can use also the ML4507 or ML4505</b></p> |                                      |    | <p><b>ML4507</b><br/>Cable and connector, <b>Snap</b>.<br/><b>ML4505</b><br/>Cable and connector, <b>Grabber</b></p> |
| <p><b>ML4866</b><br/>DRÄGER SIEMENS, Sicard 440/740 For stress Test 15 Male Pins. Cable 72829. Cardiocontrol CC1201. <b>Can be used with ML4877, ML2488 and ML4419</b> ECG cable Cardio Control Perfect MD.<br/><b>ML4877</b><br/>Cabo Précordial 10 Terminais <b>Snap</b></p>                    |  <p><b>ML4866</b> <b>ML4877</b></p> |    | <p><b>ML2488</b><br/>Leadwires <b>Banana</b><br/><br/><b>ML4419</b><br/>Leadwires <b>Clip, Snap</b>.</p>             |







|  |   |  |  |
|--|---|--|--|
| <p><b>ML7472</b><br/>Welch Allyn CP100/CP200 ECG Patient Cable 10 Lead, Resting ECG, AHA. CP-100 / CP-200 Set for Welch Allyn CP 100 / CP 200<br/><b>Patient Cable Set of 10 leads</b></p>   |                                        |   | <p><b>ML7473</b><br/>Welch Allyn Patient Lead Cable for CP100 / CP200 ECGs 400294</p>  |
| <p><b>ML7802</b><br/>Siemens / Hellige EKG cable ,10-lead ,AHA, DB15M ,need EKG leadwire ML7804</p>  |                                        |   | <p><b>ML7804</b><br/>Siemens / Hellige 412680-001, Multi-Link Leadwire set , Banana ending ,AHA</p>  |
| <p><b>ML2009</b><br/>Cabo PC-109, Kanz EKG cable , IEC , DB 15PIN &gt;DIN 3.0mm plug</p>   |                                        | <p>PC-109, Kanz EKG cable, IEC , DB 15PIN &gt;DIN 3.0mm plug.<br/>PC-109, Kanz EKG cable , IEC , DB 15PIN &gt;Banana 4.0, Kanz EKG cable , AHA , DB 15PIN &gt;Banana 4.0</p>   |  |
| <p><b>LM4488</b><br/>Cabo p/ECG Cardioline/Bosch EKG 601. ( 71251 )<br/>MONOBLOC 10 VOIES<br/><b>BANANE.</b> 10-Lead EKG Cable Cardioline Â 10-Lead Shielded EKG Cable (non-detachable leads, 4mm safety Banana), AHA American Standard<br/>This EKG cable is BRAND NEW and works for the following machines:<br/>Cardioline Delta 3 , Cardioline Delta 30 , Cardioline Delta 60 , Cardioline Delta 1 plus , Cardioline Delta 3 plus ,Cardioline Delta 60 plus</p> |                                       |    | <p><b>LM4433</b><br/>PC-104, Kanz EKG cable , IEC , DDK-16PIN &gt;DIN 3.0,</p>   |
| <p><b>LM4044</b><br/>Cabo ECG MD05F<br/>Electrocardiograph Cardio Care 2000 and brand BIONET is a monobloc cable 10 lead (banana type) for the models CardioCare and cardiotech from BIONET. Connector is SUB-D 15 (15 pins) ref : 71110. CARDIOVIT (AT-1 / AT-101/ AT-102 / AT-2 / AT-2 PLUS / AT-4 / AT-110 / AT-104 / PT-160 / Schiller AT1) Dr.LEE-310A, Dr.LEE-310A</p>   |                                      |   | <p>Philips(HP): 1772A,M3703C , M2462A. Pagewriter XLI M1700A. Dr. Lee : All HEWLETT PACKARD pagewriter 300 pi EKG. M3703A/C</p>            |
| <p><b>MK4486</b><br/>Cable for ECG PageWriter 100 Philips Mod. M1772A, 300PI<br/>Cabo ECG com 10 pontas ref. M3703C Philips. M1770A<br/>PageWriter 200i, 1772A, M3703C, M2462A, M1771A PageWriter 200; M1772A PageWriter100; PageWriter XLI M1700A. Dr.Lee: All. Cardiorapid K33-M, Cardiorapid K33M</p>   |                                     |   | <p><b>MK4486</b><br/>Cabo ECG com 10 pontas ref. M3703C Philips. PageWriter 300PI</p>  |
| <p><b>ML5844</b><br/>PageWriter 300 PI Mod. M1770A. •M3703C/M3703A PageWriter 300PI.<br/>M3703C/M3703A: Esaote Replacement patient cable<br/><b>HP EKG CABLE and leadwires , IEC, Grabber</b></p>  |  <p><b>ML5844</b> <b>ML2001</b></p> | <p><b>ML2001</b><br/>PageWriter 300 PI Mod. M1770A. •M3703C/M3703A PageWriter 300PI.<br/>M3703C/M3703A: Esaote Replacement patient cable. Contec 300G. Contec 600G<br/><b>HP EKG CABLE and leadwires , IEC, Banana 4.0</b></p> | <p><b>ML5845</b><br/>PageWriter IEC, <b>SNAP</b></p>  |















|   |  |  |   |
|---|--|--|---|
| <p><b>RC444130</b><br/>Nihon Kohden BJ-902D (K081) (BJ-902D) BJ 902D, ECG 9010, 9020, 9022, 9110, 9130, 9132, ECG 9620, 1250, 1350, PEA-1110. Plug K081<br/>ECG Nihon Fukuda Cardiofax GEM<br/><b>Original Nihon Kohden Cable</b></p>   |     | <p>Connection cable for 10 electrodes. Color Coded 4mm diam plug with spring contact (IEC). Cable length: limb ca. 3.4m. chest ca 2,9m.<br/><b>Original Nihon Kohden Part Number: BJ-902D / K081.</b></p>  |   |
| <p><b>RC44488</b><br/>Nihon Kohden 10-lead EKG cable, AHA, NK DB15M &gt;banana 4.0, No Resistor.<br/><b>Cardioline:</b> E1.E1R. E3, Delta 1, 300, 360, C3</p>   |     | <p><b>Nihon Kohden</b> ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, <b>DongJiang</b> ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120. BIOCARE ECG-101, ECG-101G, ECG-300, ECG-300G, ECG-1200, ECG-9803</p>  |   |
| <p><b>RC88466</b><br/>BJ- 901D, Nihon Kohden 10-lead EKG cable ,IEC ,NK DB15M &gt;DIN3.0 ,No resistor<br/><b>RC88464</b><br/>Nihon Kohden 10-lead EKG cable ,IEC ,NK DB15M &gt;DIN3.0 , <b>4.7KΩ resistor</b></p>   |     | <p><b>Nihon Kohden</b> ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, <b>DongJiang</b> ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120. BIOCARE ECG-101, ECG-101G, ECG-300, ECG-300G, ECG-1200, ECG-9803</p>  | <p><b>Cardioline;</b> E1.E1R. E3, Delta 1, 300, 360, C3</p>                                     |
| <p><b>ML6642</b><br/>FUKUDA ALPHA 600AX Cable ECG KP-310, FUKUDA ME CX 600 AX Fukuda KP-310, 501B-III FUKUDA DENSHI 501B III single ch<br/><b>Old Model</b><br/><b>Round 12pin&gt;Banana 4.0</b></p>  |    |   |             |
| <p><b>RC44413</b><br/>Nihon Kohden 10-lead EKG cable ,AHA ,NK DB15M &gt;banana 4.0 , <b>4.7K resistor</b><br/>Nihon Kohden BJ-902D, K081, Connection cable for 10 electrodes · colour, coded 4 mm diam plug with spring contact (IEC) · cable length: limb ca. 3.4 m, chest ca. 2.9 m. Cardiofax. ECG 3 pistes NIHON KOHDEN 1150. I.P. 72258. BJ 902D, BJ902D,</p>  |  | <p>Nihon Kohden: Cardiofax 6151 /6353, Cardiofax EKG-8270/8350, 8370/8420 /8423/8450/8453/8820/8823/9130, Cardiofax-O. Fukuda Cardisunny alpha 1000, Cardisunny alpha 6600A. ECG-9620/1250/1350. Fukuda KP-313. Nihon Kohden Cardiofax S. Fukuda Alpha 600AX, KP-310 (K310). Plug K081 ECG Nihon Fukuda Cardiofax GEM</p>  | <p>ECG Siemens (Megacart). ECG AT-3 Shiller ref. 3404.7501.0 0. Cabo ECG 9020K Nihon Kohden</p> |
| <p><b>ML4092</b><br/>Fukuda Cardisunny Alpha 1000 (Cardisunny Alfa 1000). Fukuda 501B Fukuda KP-313 (Fukuda KP313) Nihon Kohden one-piece EKG cable with 10-lead NK DB15M &gt;banana 4.0, 4.7K resistor. Cardiofax 6151/6353, Cardiofax EKG-8270/ 8350/ 8420/ 8370/ 8423/ 8450/ 8453/ 8820/ 8823/ 9130 Cardiofax-Q, FUKUDA ME: (Cardisunny Alfa 1000) Cardisunny alpha 1000, Cardisunny alpha 6600A/ 6600 AX. Burdick: EK-10, Elite, Elite II, E-350, E350i, E550, E560, (1st Generation Eclipse LE/Lei, with serial number starting with 1051). Burdick#007705, Stress cable with detachable leads 007706.</p> |   | <p>Cardisunny alpha 6600A/ 6600 AX Burdick: EK-10, Elite, Elite II, E-350, E350i, E550, E560, (1st Generation Eclipse LE/Lei, with serial number starting with 1051). Burdick#007705, Stress cable with detachable leads 007706</p>  <p><b>NK DB15M &gt; DIN 3.0, 4.7K resistor.</b><br/>Connector: DB 15pin male, 3.0 Din type needle on electrode end<br/><b>4.7 KL Resistance</b><br/><b>See RC88464</b></p> | <p>Fukuda Cardisunny 501-B III. Fukuda Cardisunny 501B III</p>                                  |
| <p><b>ML5402</b><br/>FUKUDA ME: Cardisunny Alpha 1000 (Cardisunny Alfa 1000). Nihon Kohden: Cardiofax 6151 /6353, Cardiofax EKG-8270/8350, 8370, 8420 / 8423/ 8450/ 8453/ 8820/ 8823/ 9130, Cardiofax-O Fukuda Cardisunny Alpha 1000, Cardisunny Alpha 6600A. (Nihon Kohden ECGLAB Treadmill test cable, 10-Lead EKG cable IEC, DB15M-&gt;Snap, 4.7K resistor.</p>  |   |   | <p><b>ML5404</b><br/>10-Lead EKG cable ,IEC, NK DB15M-&gt;Grabber (CLIPS), 4.7K resistor</p>    |




|   |   |  |   |
|---|---|--|---|
| <p><b>ML1018</b><br/>Fukuda Cardisuny Alpha 1000<br/>Fukuda ME KP-500D EKG<br/>Cable , IEC, DB15M&gt;<b>Grabber</b><br/><b>See that the position of the<br/>connector is reversed of the<br/>ML1010</b><br/><b>We Have this cable with banana too</b></p>               |    |   |    |
| <p><b>ML1010</b><br/>Fukuda ME KP-500D EKG cable<br/>,IEC, DB15M&gt; <b>Banana 4.0</b><br/>Fukuda ME Cardisuny KP-500<br/>Cardisuny Fukuda KP-313<br/>Fukuda Cardisuny Alpha 1000<br/>Fukuda CX1000, Fukuda CX 1000</p>   |    |   | <p>Fukuda<br/>CX1000,<br/>Fukuda ME<br/>Cardisuny<br/>C-110</p>   |
| <p><b>ML1011</b><br/>Fukuda ME KP-500D EKG Cable<br/>,IEC, DB15M&gt;<b>Grabber</b> ,<br/>Fukuda ME Cardisuny KP-500<br/>Cardisuny Fukuda KP-313<br/>Fukuda Cardisuny Alpha 1000 Fukuda<br/>CX1000, Fukuda CX 1000</p>   |   | <p><b>ML1098</b><br/>Fukuda ME KP-500D EKG<br/>Cable ,IEC, DB15M&gt;<b>SNAP</b><br/>20KΩ resistor<br/>Fukuda ME Cardisuny KP-500<br/>Cardisuny Fukuda KP-313<br/>Fukuda Cardisuny Alpha 1000</p>   |    |
| <p><b>ML1197</b><br/>Fukuda Denshi <b>5-lead Graber</b> 15pin<br/>ECG cable for DS-8000 Series<br/>Monitors.<br/>Fukuda Dynascope 8000 serie<br/>Denshi DS-8000 series: DS-8100,<br/>DS-8500, DS-8007, DS-8200, DS-8400</p>   |  | <p><b>ML1196</b><br/>Fukuda Denshi <b>5-lead Snap</b><br/>15pin ECG cable for DS 8000<br/>Series Monitors. Fukuda<br/>Dynascope 8000 serie. Denshi<br/>DS-8000 series: DS-8100, DS-8500,<br/>DS-8007, DS-8200, DS-8400</p>                               |   |
| <p><b>ML1977</b><br/>Fukuda multi-link ECG trunk cable,<br/>10-lead, IEC, Round 12pin, 20K resistor.<br/>Fukuda beside <b>monitor DS-7300</b><br/>series with 12pin connectors.</p>  |  | <p><b>ML1978</b><br/>ECG leadwire, 5-Lead ,Snap<br/>,IEC, 40IN &amp; 2pin, connect trunk cable<br/><b>N,C,R,F,L</b></p>  <p><b>We have other configurations</b></p> | <p><b>ML1979</b><br/>ECG leadwire<br/>set,5-<br/>LeadSnap, IEC,<br/>40IN connect<br/>trunk cable<br/><b>C2-C6</b></p>  |
| <p><b>ML1919</b><br/>Fukuda ME FX-031 EKG cable ,IEC,<br/>Round 11pin &gt; <b>Banana 4.0</b> , No<br/>resistor<br/>Fukuda FX-031, Fukuda FX031</p>  |  | <p><b>ML1918</b><br/>Fukuda ME FX-031 EKG cable<br/>, IEC, Round 11pin &gt;<br/><b>DIN 3.0mm</b>, No resistor<br/>Fukuda FX-031 Fukuda FX031</p>   |   |
| <p><b>ML5799</b><br/>Fukuda Denshi<br/>One piece ECG cable with leadwire,<br/>5-lead, IEC, Snap, 12pin. Fukuda<br/>Denshi DS5100E, DS5300W,<br/>DS5400, DS7101/L, DS7141,<br/>DS7001. FUKUDA Compatible One<br/>piece ECG cable with leadwires</p>                      |  |   | <p><b>ML5798</b><br/>Fukuda Denshi,<br/>One piece ECG<br/>cable with<br/>leadwire, 3-<br/>lead, AHA,<br/>Grabber, 12pin.</p>  |
| <p><b>ML5818</b><br/>Fukuda Denshi, Fukuda 3-lead trunk<br/>cable, With simens 3-lead yoke, round 12pin<br/>-2PIN. Fukuda Denshi DS5100E, DS5300W,<br/>DS5400, DS7101/L, DS7141, DS7001.</p>  |  |   | <p><b>ML5772</b><br/>Fukuda Denshi,<br/>3-lead IEC<br/>Snap ECG<br/>leadwire .IEC,<br/>snaps.</p>   |
| <p><b>ML8502</b><br/>YY-1063, ANEC EKG cable ,<br/>IEC , Round-15pin &gt;<b>DIN3.0</b><br/><b>ML8503</b><br/>NEC EKG cable ,AHA , Round-<br/>15pin &gt;<b>Banana4.0</b></p>   |  | <p>NEC: K-3321, 3330, 3370<br/>NEC K-3321, NEC 3330,<br/>NEC 3370</p>  |   |





















|   |  |  |  |
|---|--|--|--|
| <p><b>MK4487</b></p> <p>Cable for ECG Philips Mod. M1772A (Cable M3703C). Compatible with the following ECG machines: EDAN: Smart SE-1, SE-3, Esaote P80, Esaote P8000, GUIMA: CARDIOGIMA 1M, 3M, 6M, 12, 12M, H.P.: Pagewriter 1513, 154B/C, 1514D, 1515A/B, 1517A., KONTRON: EKG P80, Macquarie: CT100, CT200., NIHON KOHDEN: Cardiofax 651, 6581, 8110., PROGETTI: all models, SHILLER: AT/1, AT/2, AT/3, AT/4, AT/5, AT/6, AT/10, CS6, CS12, CS100, CV6, P80., Bionet: Cardiocare, Cardiotouch, Cardioline: E1, E1R, E3, Delta 300, Delta 360 C3. CARDIETTE: (Elettronica Trentina): Autoruler 12/1, Autoruler 12/3, Excel 103, Excel 106, Astron Comp 1, Astron Comp 3, Astron Comp 6, Start 100 P, Start 100 H, Start 200, Start 200HV, AR600, AR 600 ADV, AR1200, AR1200, ADV. Cardioraoid Camina: K111, K112, K33, K36, K300 M, K300 P, K600 P, K600 D, K130 B, K 130 D. Cable F6726R.<br/>Dr. LEE ECG-310A, Dr. LEE 310A<br/>Dr. LEE-310A</p>  |           |    | <p>Esaote P 80,<br/>Esaote P 8000<br/>Dr. LEE<br/>ECG-310A<br/>Dr. LEE 310A<br/>Dr. LEE-310A</p>   |
| <p><b>RC44413,</b></p> <p>Nihon Kohden 10-lead EKG cable ,AHA ,NK DB15M &gt;<b>banana 4.0</b> ,<br/><b>4.7K resistor</b></p> <p>Nihon Kohden BJ-902D, K081, Connection cable for 10 electrodes - colour, coded 4 mm diam plug with spring contact (IEC) - cable length: limb ca. 3.4 m, chest ca. 2.9 m. Cardiofax. ECG 3 pistes NIHON KOHDEN 1150. Integral process 72258. BJ 902D, BJ902D<br/>ECG Nihon Fukuda Cardiofax GEM</p>  |     | <p>Nihon Kohden:<br/>Cardiofax 6151 /6353,<br/>Cardiofax EKG-8270/8350,<br/>8370/8420<br/>/8423/8450/8453/8820/8823/91<br/>30, Cardiofax-O. Fukuda<br/>Cardisuny alpha 1000,<br/>Cardisuny alpha 6600A. ECG-<br/>9620/1250/1350. Fukuda KP-<br/>313. Nihon Kohden Cardiofax<br/>S. Fukuda Alpha 600AX<br/>KP-310 (K310). Plug K081</p> | <p>ECG<br/>Siemens<br/>(Megacart).<br/>ECG AT-3<br/>Shiller ref.<br/>3404.7501.0<br/>0. Cabo<br/>ECG 9020K<br/>Nihon<br/>Kohden</p>  |
| <p><b>MK4485</b></p> <p>ESAOTE, COD. F6725R-F6725RP-F6200R-F6200RP-F6209-F6209P. 'LOW MOISE' 10 LEADS PATIENT CABLES FOR THE FOLLOWING ECG DEVICES:<br/>ESAOTE: model: F6725R - 10 leads patient cable with 4mm plug, model: F6725RP - 10 leads patient cable with clamps, model: F6200R - 10 leads patient cable with 4mm plug, model: F6200RP - 10 leads patient cable with clamps, model: F6209 - 10 leads patient cable - round connector with 4mm plug, model: F6209P - 10 leads patient cable round connector with clamps. Compatible with the following ECG machines: EDAN: Smart SE-1, SE-3, ESAOTE: P80, P8000, GUIMA: CARDIOGIMA 1M, 3M, 6M, 12, 12M, H.P.: Pagewriter 1513, 154B/C, 1514D, 1515A/B, 1517A., KONTRON: EKG P80, Macquarie: CT100, CT200., NIHON KOHDEN: Cardiofax 651, 6581, 8110., PROGETTI: all models, SHILLER: AT/1, AT/2, AT/3, AT/4, AT/5, AT/6, AT/10, CS6, CS12, CS100, CV6, P80., Bionet: Cardiocare, Cardiotouch, Cardioline: E1, E1R, E3, Delta 300, Delta 360 C3. SHILLER AT1, AT2, AT3, AT4, AT5, AT6, AT10</p> |    |  <p>Plug K081</p>   | <p>CARDIETTE:<br/>(Elettronica Trentina):<br/>Autoruler 12/1,<br/>Autoruler 12/3,<br/>Excel 103,<br/>Excel 106,<br/>Astron Comp 1,<br/>Astron Comp 3,<br/>Astron Comp 6,<br/>Start 100 P,<br/>Start 100 H,<br/>Start 200, Start<br/>200HV, AR600,<br/>AR 600 ADV,<br/>AR1200,<br/>AR1200, ADV.<br/>Cardioraoid<br/>Camina: K111,<br/>K112, K33,<br/>K36, K300 M,<br/>K300 P, K600<br/>P, K600 D,<br/>K130 B, K 130<br/>D</p> |




|   |   |   |   |
|---|---|---|---|
| <p><b>ML5609</b><br/>Bruker / Physiogard 910<br/>Bruker One-piece ECG cable, 4-<br/>Lead , Snap , AHA,12ft, Round<br/>12pin<br/>Bruker One-piece ECG cable, 4-<br/>Lead , Snap , IEC 12ft, Round 12pin</p>  |    |            |   |
| <p><b>ML5618</b><br/>Schiller 2.400209<br/>3-lead ECG cable with 3-lead<br/>IEC grabber leadwires.</p>  |    | <p>Part n. 2.400209, Schiller<br/>3-lead ECG cable with 3-<br/>lead IEC grabber leadwires</p> |   |
| <p><b>ML9112</b><br/>HP/Philips, M1663A, HP 10-lead<br/>Trunk cable , HP <b>12pins .2M</b><br/>Philips. 78352C, 78354C, 78834C, Legacy<br/>HP patient MonitorsHPA5-90P, HPB5-90P,<br/>HP A5-90P, HP B5-90P,<br/>862472 to 862475 C3 Monito862478,<br/>862479 C3 Monitor,863064, 863065,<br/>863066, 863068<br/>SureSigns VM6, VM8 ECG Modules ,<br/>M1001A/B, M1002A/B, M3000A, M3001A,<br/>M3002A MMS X2 , M8102A MP2 Monitor ,<br/>M8105A MP5 Monitor<br/>M4735A HeartStart XL defibrillator monitor ,<br/>M3535A/M3536A HeartStart MRx</p> |    |            | <p><b>ML9114</b><br/>HP/Philips,<br/>M1949A, HP<br/>10-lead<br/>Trunk cable ,<br/>HP<br/><b>12pins .27M</b><br/>HPA5-90P,<br/>HPB5-90P,<br/>HP A5-90P,<br/>HP B5-90P,</p> |
| <p><b>ML9113</b><br/>HP/Philips M1971A+M1978A, HP<br/>10-lead ECG leadwire, Grabber,<br/>AHA. HPA5-90P, HPB5-90P, HP<br/>A5-90P, HP B5-90P, Philips. 78352C,<br/>78354C, 78834C, Legacy<br/>862472 to 862475 C3 Monito862478,<br/>862479 C3 Monitor 863064, 863065,<br/>863066, 863068 SureSigns VM6, VM8 ECG<br/>Modules , M1001A/B, M1002A/B, M3000A,<br/>M3001A, M3002A MMS X2 , M8102A MP2<br/>Monitor , M8105A MP5 Monitor, M4735A<br/>HeartStart XL defibrillator monitor ,<br/>M3535A/M3536A HeartStart MRx</p>                        |  |          | <p><b>ML9123,</b><br/>HP/Philips<br/>M1645A+M1<br/>604A,<br/>HP 10-lead<br/>ECG<br/>leadwire.<br/>Snap, IEC<br/>HPA5-90P,<br/>HPB5-90P,<br/>HP A5-90P,<br/>HP B5-90P,</p> |
| <p><b>ML8707</b><br/>CU Medical 10-lead ECG cable<br/>with leadwire , IEC Round<br/>12pin-&gt; Snap</p>   |  |          |   |
| <p><b>ML5401</b><br/>Nihon Kohden 10-Lead EKG cable,<br/>IEC, NK DB15M-&gt;Grabber, No<br/>resistor. ECG-9132, DongJiang<br/>ECG-11A, ECG-11B, ECG-11D,<br/>ECG-32A, (II), ECG-120. Biocare<br/>ECG-101, ECG-101G, ECG-300,<br/>ECG-300G, ECG-1200, ECG-9803,<br/>Cardioline; E1.E1R. E3, Delta 1,<br/>300, 360, C3. Carewell</p>   |  |          | <p><b>ML5403</b><br/>Nihon<br/>Kohden<br/>ECGLAB<br/>Traadmill<br/>test cable ,<br/>IEC,<br/>DB15M-<br/>&gt;Snap , No<br/>resistor</p>                                    |
| <p><b>ML8788</b><br/>Nihon Kohden BR-912D, Nihon<br/>Kohden ECG-9320 10-lead EKG<br/>cable, IEC. Use with patient<br/>input-box JC-901D for resting ECG</p>   |  |          | <p>Use with<br/>patient input-<br/>box JC-901D<br/>for resting<br/>ECG</p>  |





|  |   |   |   |
|--|---|---|---|
| <p><b>ML7542</b><br/>Charmcare, spo2 adapter cable , Redel 6pin with 1key -&gt; DB9, 2.5M<br/><b>Use with Charmcare Model PRIZM3</b></p>   |    |    | <p>Use with Charmcare Model PRIZM3</p>  |
| <p><b>ML2698</b><br/>GE Cable 2016560-002, ECG. CAM 14, Coiled, Patient, Round Connector. Length <b>2,70/4,70m</b>. 2000 Treadmill CAM 14 Coiled cable, GE/Marquette Case, Case 8000, Mac 5500, 5500ST, 5000.</p>                          |   |    | <p><b>ML2694</b><br/>Cable with <b>4m/6m</b><br/><b>ML2696</b><br/>Cable with <b>6m/8m</b></p>                              |
| <p><b>ML5834</b><br/>GE MAX 1, 700044-101, GE AM4/AM5 Coiled interface cable, 15ft for MAX 1. Use with GE MAX 1</p>  |    |    | <p><b>ML5833</b><br/>GE 700044-204 GE AM4/AM5 Coiled Interface Cable for GE MAX-Personal, MAC 8 Systems with AM4 or AM5</p> |
| <p><b>RC554910</b><br/>ShangHai Kohden EKG Cable , DB15M&gt;<b>DIN 3.0</b></p>   |    | <p><b>ShaiHai Kohden</b> : ECG-6511, ECG-6151, 8110P/K, 6353,</p>    | <p><b>ShangHai Kohden</b></p>   |
| <p><b>ML1159, RC55422B</b><br/>ShangHai Kohden EKG Cable, DB15M&gt;<b>Banana 4.0</b>. Cardioline AR1200ADV. CARDIOLINE AR600 AR1200 AR2100</p>   |   | <p><b>ShaiHai Kohden</b> : ECG-6511, ECG-6151, 8110P/K, 6353,</p>   | <p><b>ShangHai Kohden</b></p>   |
| <p><b>ML6622</b><br/>MAC-400/MAC-500/MAC-1200 One-piece EKG CABLE and leadwires, IEC, AHA <b>Banana</b><br/>Hellige: MicroSmart MAC 400, MAC 500, MAC1100, MAC 1200/1200 ST. Cabo para Eletrocardiógrafo GE Modelo MAC 800. MAC 1200ST</p> |  | <p>GE ECG Machine Parts 10 Lead Cable LDWR IEC 2104726-001 Medical Device.</p>   | <p>Hellige: MicroSmart MAC400, MAC500, MAC800, MAC1200 MAC1200ST</p>  |
| <p><b>ML6634</b><br/>MAC-400/MAC-500/MAC-1200 One-piece EKG CABLE and leadwires, IEC, AHA <b>SNAP</b><br/>Hellige: MicroSmart MAC 400, MAC 500, MAC1100, MAC 1200/1200 ST. Cabo para Eletrocardiógrafo GE Modelo MAC 400. MAC 1200ST</p>   |  |    | <p><b>ML6633</b><br/>MAC-400, MAC-500 MAC-1200 One-piece EKG CABLE and leadwires, IEC, AHA <b>Grabber</b></p>               |
| <p><b>ML5156</b><br/>Multi-Link <b>Cable</b>, IEC, Mac-500/ MAC-1200ST , GE Marquette MAC500, MAC1200. Mac 500/ MAC 1200 ST MAC 1200ST.<br/><b>Use with ML5157 or ML5159</b></p>   |  |    | <p>MAC1200 MAC1200ST</p>  |
| <p><b>ML5157</b><br/>Multi-LiNK <b>Leadwire set</b>, <b>Banana</b> ending, IEC, Mac1200ST, Cable 412680-001. Mac 500/ MAC 1200 ST MAC 1200ST.<br/><b>Use with ML5156</b></p>   |  | <p><b>ML5159</b><br/>Multi-LiNK <b>Leadwire set</b> , <b>Grabber</b> ending, IEC, Mac1200ST Mac 500/ MAC 1200 ST MAC 1200ST.<br/><b>Use with ML5156</b></p>  |                                        |




|   |  |   |  |
|---|--|---|--|
| <p><b>ML7748</b><br/>Multi-LiNK Leadwire set ,<br/><b>Banana ending,IEC</b></p> <p><b>ML7748A</b><br/>Multi-LiNK Leadwire set ,<br/><b>Snap ending,IEC</b><br/>2000 Treadmill</p>   |    |    | <p><b>ML7748B</b><br/>Multi-LiNK<br/>Leadwire<br/>set. 10<br/>900178-003<br/><b>Grabber<br/>ending,IEC</b></p> |
| <p><b>ML2698</b><br/>GE Cable 2016560-002, ECG.<br/>CAM 14, Coiled, Patient, Round Connector.<br/>Length <b>2,70/4,70m</b>. 2000 Treadmill<br/>CAM 14 Coiled cable, GE/Marquette Case,<br/>Case 8000, Mac 5500, 5500ST, 5000.</p> |    |    | <p><b>ML2694</b><br/>Cable with<br/><b>4m/6m</b><br/><b>ML2696</b><br/>Cable with<br/><b>6m/8m</b></p>         |
| <p><b>ML2699</b><br/>GE Marquette CAM 14 HD<br/>Module<br/>GE Module 900995-002, CAM 14.<br/>F3MF1916P</p>  |    |    | <p>2000<br/>Treadmill</p>  |
| <p><b>ML5648</b><br/>Hellige EK51, 10-lead EKG<br/>cable with <b>banana 4.0mm</b>.<br/>Hellige EK-51, Hellige EK 51<br/>Contec 300G, Contec 600G</p>  |    | <p>Hellige EK-51, Hellige EK 51<br/>10-lead EKG cable ,IEC,<br/>AHA, NK DB15M &gt;banana<br/>4.0</p>  | <p>Cabo ECG<br/>para<br/>Hellige EK51.<br/>Cabo ECG<br/>para<br/>Hellige EK 51</p>                             |
| <p><b>RL8893</b><br/>BC-753V K351A Nihon Kohden one<br/>piece ECG cable with leadwire3,-<br/>lead, IEC, Grabber, 11pin.<br/>Nihon Kohden Mod. TEC-7731K<br/>Nihon Kohden TEC-7631K</p>  |   | <p>Cardiolife Defibrillators TEC-<br/>5500, TEC-7600 , TEC-7700 ,<br/>BSM-7102 , OMP-7201, OMP-<br/>7203 , TEC-5200A ,<br/>(S/N&lt;20421) <b>BC-753V K351A</b><br/>Nihon Kohden TEC-7731K</p> |                            |
| <p><b>RL8893A</b><br/>BC-763V K350A, Nihon Kohden<br/>one piece ECG cable with<br/>leadwire3,-lead, IEC, Snap, 11pin.<br/>Nihon Kohden Mod. TEC-7731K<br/>Nihon Kohden TEC-7631K</p>  |  | <p>Cardiolife Defibrillators TEC-<br/>5500, TEC-7600 , TEC-7700 ,<br/>BSM-7102 , OMP-7201, OMP-<br/>7203 , TEC-5200A ,<br/>(S/N&lt;20421) <b>BC-763V K350A</b><br/>Nihon Kohden TEC-7731K</p> |                           |
| <p><b>RL8895</b><br/>Nihon Kohden one piece ECG cable<br/>with leadwire, 5-lead, AHA, Snap<br/>11pin.<br/>Nihon Kohden Mod. TEC-7731K<br/>Nihon Kohden Mod. TEC-7631K</p>   |  | <p>Cardiolife Defibrillators TEC-<br/>5500, TEC-7600 , TEC-7700 ,<br/>BSM-7102 , OMP-7201, OMP-<br/>7203 , TEC-5200A ,<br/>(S/N&lt;20421)<br/>Nihon Kohden TEC-7731K</p>                      |                           |
| <p><b>RL4493</b><br/>Nihon Kohden trunk cable,<br/>11pin. Nihon Kohden Mod.<br/>TEC7731K<br/>Nihon Kohden TEC-7631K</p>   |  | <p>Cardiolife Defibrillators<br/>TEC-5500, TEC-7600 ,<br/>TEC-7700 , BSM-7102 ,<br/>OMP-7201, OMP-7203 ,<br/>TEC-5200A , (S/N&lt;20421)</p>   |                           |
| <p><b>RL15912</b><br/>one piece ECG cable with<br/>leadwire3,-lead,AHA,Snap,Metal<br/>Round 8 pin</p>   |  | <p>Nihon Kohden OEC-6102A<br/>(s/n&lt;20156) , OEC-6105A<br/>(s/n&lt;20226) , OEC-6120A ,<br/>OEC-6120J , OEC-6120K</p>   | <p>All Nihon<br/>Kohden<br/>Compatible<br/>8pin will use<br/>Metal 8PIN</p>                                    |
| <p><b>RL15918</b><br/>One piece ECG cable with<br/>leadwire5,-lead,IEC,Grabber<br/>,Metal Round 8 pin. All Nihon<br/>Kohden Compatible 8pin will use<br/>Metal 8PIN</p>   |  | <p>Nihon Kohden OEC-6102A<br/>(s/n&lt;20156) , OEC-6105A<br/>(s/n&lt;20226) , OEC-6120A ,<br/>OEC-6120J , OEC-6120K</p>   | <p>All Nihon<br/>Kohden<br/>Compatible<br/>8pin will use<br/>Metal<br/>8PIN</p>                                |




|   |   |   |   |
|---|---|---|---|
| <p><b>ML2348</b><br/> Module Pagewriter Trim III<br/> <b>Patient Interface Module.</b><br/> Philips Pagewriter Trim Series Patient Module 989803130101, 989803130091 For Philips M1700A, M1701A, M1716B, M1702A Trim II, Trim 1, 2, 3.<br/> <b>Without 10 Leads Set ECG Cable</b></p> <p><b>ML2348A</b><br/> <b>Complete Lead Wire Set and Data Cable</b></p>   |    <p><b>ML2348</b>      <b>ML2341</b></p>   | <p><b>ML2341</b><br/> USB Data wire for HP trim ECG box. Data wire for HP trim ECG acquisition box; cable length: 1.8 meters; Compatible with: Pagewriter Touch (860284) – Pagewriter Trim I , II , III (860286-860287-860288-860290) Use with HP trim leadwire set. Phillips Trim USB Cable 453564034571</p>  |   |
| <p><b>ML2331</b><br/> Philips 10 Leads Set Ecg Cable For Philips M1700A, M1701A, M1716B, M1702A Trim II, Trim 1, 2, 3 Philips PageWriter Trim III, Trim II, Trim I , M1702A, M1713B Cable 989803130101 Leadwires 989803130101</p>   |      | <p>Complete Lead 10 color, IEC, Banana 4.0 Set for Philips Trim, Trim I, Trim II (989803130101), Trim III, PageWriter TC30, PageWriter TC50, PageWriter XLI M1700A, PageWriter XLE 1702A, Philips M1701A, M1702A; Trim II, Trim 1, 2, 3 Philips PageWriter Trim III, Trim II, Trim I</p>  |    |
| <p><b>ML2332</b><br/> Cable M1716B, HP M1716B EKG Leadwire set, IEC, 2PIN &gt; <b>Grabber.</b><br/> Philips PageWriter Trim III, Trim II, Trim I. Complete lead sets, 10 color marked rings, IEC, <b>Grabber</b> Philips Trim Touch Leadwires 989803129161</p>  |     | <p>Philips 10 Leads Set Ecg Cable For Philips M1700A, M1701A, M1702A Trim II, Trim 1, 2, 3 Philips PageWriter Trim III, Trim II, Trim I Hp ECG Leadwire set ,IEC , 2pin -&gt; Grabber</p>   | <p><b>ML2333.</b><br/> EKG leadwire set, IEC, 2PIN &gt; <b>Snap</b></p>   |
| <p><b>ML8490</b><br/> Philips PageWriter TC30, TC50 TC70, 989803151641, Philips 10-lead leadwires, IEC, Plastic 14pin&gt;Banana 4.0, <b>Limb lead 1m.</b> Philips 989803151641 Pagewriter TC70 <b>TC Series Complete Leadwire Set</b></p>   |      |     | <p><b>ML2333.</b><br/> EKG leadwire set, IEC, 2PIN &gt; <b>Snap</b></p>    |
| <p><b>ML4155</b><br/> Cables for HP. M1510A 71780S EMBASE 3 VOIES PHILIPS-AGILENT. H.P. 12 CONTACTS MALES IEC</p>   |     | <p><b>ML6678</b><br/> Cables for HP. M1510A EMBASE 3 VOIES PHILIPS-AGILENT. H.P. 12 CONTACTS MALES IEC</p>  |    |
| <p><b>RL8421</b><br/> Zoll 579-A-3ST. MONOBLOC <b>3 LEADS SNAP ECG AHA,</b> SPACELABS, SCOPE 90701-90303-90305 –90308 = PC EXPRESS90309 = PC SCOUT90385 = UCW 6 MALE PINS PROTECTED. Contec CMS6000</p>    <p>Criticare Ref<sup>®</sup> 1073/IS 83305B002. RL8421B,</p> | <p><b>RL8421A</b> , One piece ECG cable with 3-ld , TPU, <b>ECG IEC. Clip,</b> 12ft &amp; 180° straight connector , Use with Instrument : BCI :3101,3404,6100,6200,9100,9200, Advisor , Autocorr Colin Medical: BP308, BP408, CBM3000 CSI : 504US, 602-4usp, poet 2, Scholar GE Medical Corometrics 556, Critikon 8100T, Critikon Dinamap 8700 Series, Critikon Dinamap 9700 series. Patient Monitor Contec CMS-6000 MEDTRONICS PHYSIO-CONTROL : Lifepak 6, VSM-1ESF(w/resistor ); NELLCOR :N-200 C-LOCK,N-250 C-clock,N-1000,NPB-3900,NPB-4000, Nihon Kohden :BSM-8301A, BSM-8500A, BSM-8502A ,OEC-6102A/6105A/7102A,8108A, TEC-5200a, TEC-7100A/7200A/7300A;Pace tech Minipack 300,Minipack 3000/3100, Vitalmax 800 Plus , Vitalmax 2200, Vitalmax 4000, Vitalmax4000 Modular, Vitalmax 4100,Vitalmax 4100CL; Philips Medical Systems : M3923A A1 Compact Monitor, M3925 A1 Compact Monitor, M3927 A3 Compact Monitor, M3929A A3 Compact Monitor, Welch Allyn: PROPAQ 100 series, Propaq 200 Series, Propaq 200 Series, Propaq Encore. Contec CMS6000/ Contec CMS7000/CMS8000/CMS9000</p> |      | <p><b>RL8421</b><br/> 1073/S 3 Lead ECG Cable, Criticare, <b>3 Lead Snap ECG AHA</b> Cable Assembly, ECG / EKG - ECG / EKG Cables Criticare Part #: 1073/S , SKU #: FME-64357 , 3 Lead Snap ECG AHA Cable Assembly (8100/E/H)</p>  |
























|   |  |  |  |
|---|--|--|--|
| <p><b>ML5672</b><br/>Cabo ERBE Microwave ERBE<br/>12-100,12-240, UHF 69B, 1100 B<br/>Erbotherm ThermoTherapy</p>    |   | <p><b>ML1445</b><br/>Válvula 3-750 "Valvo" P/Ondas<br/>Curtas, Short-Wave<br/>Válvula TB 3/750 (RS 1026)<br/>Erbotherm ThermoTherapy</p>      | <p>Fusível H. T.<br/>500mA<br/>(80x6,2mm)<br/>ERBE</p>  <p><b>ML5524</b></p>                      |
| <p><b>ML1431</b><br/>ERBE<br/>MAGNETRÃO (Magnetron)<br/>Magnetron VALVO 7090A<br/>Microwave ERBE 12-100,12-240<br/>Erbotherm 12-240, Erbotherm 12-100<br/>Erbotherm ThermoTherapy</p>   |  <p><b>ML1431</b></p>   |  <p><b>ML6506</b></p>   | <p><b>ML6506</b><br/>Ponte<br/>rectificadora<br/>H.T.<br/>H441-1<br/>(VARO 8205)<br/>Erbotherm<br/>12-100, 12-<br/>240</p>   |
| <p><b>ML6605</b><br/>Relé NMR 2hom 400w ERBE, Erbotherm</p>    | <p><b>ML6507</b><br/>Bobine H.F. ERBE, Erbotherm ThermoTherapy</p>  | <p><b>ML4427</b><br/>Reóstato R100 27Kohm 0.06A<br/>2/84 (Diâmetro 84mm) ERBE</p>    | <p><b>ML6511</b><br/>Amperímetro<br/>400MA ondas<br/>Curtas ERBE<br/>12-240<br/>AQR85 400mA</p>  |
| <p><b>ML2181</b><br/>Cables for Defibrillator Physio-<br/>control Life pack 9, Ref. 800947-01<br/>One piece ECG cable with 3-ld,<br/>TPU, IEC, Snap, 12ft &amp; 120° angle<br/>connector, 6pin. 008-0315-01<br/>ZOLL M Series, Lifepack 9</p> |  <p><b>ML2181</b>      <b>ML2182</b></p>                           | <p><b>ML2182</b><br/>Cables for Defibrillator Physio-<br/>control Lifepack 9<br/>One piece ECG cable with 3-ld, TPU,<br/>IEC, Clip, 12ft &amp; 180° straight<br/>connector .6pin . 008-0315-01<br/>ZOLL M Series, Lifepack 9</p> |   |
| <p><b>RL6661</b><br/>MEDTRONIC<br/>PHYSIO CONTROL LIFEPAK 9<br/>Cabo desfibrilhador lifepak 9<br/><b>5 LEADS CLIP</b>, 72507</p>  |    | <p><b>ML5014</b><br/>800947-01 MEDTRONIC<br/>PHYSIO Control LIFEPAK 9<br/>Cabo desfibrilhador Lifepak 9<br/><b>3 LEADS CLIP</b>. 72504</p>   | <p>LIFEPAK 5-<br/>6-7-8-10<br/>VSM1-2-3-<br/>1440, 6<br/>MALE PINS</p>   |
| <p><b>ML6410</b><br/>ECG cable EC6410 Mindray.<br/>Original Cable.<br/>MINDRAY - Mod. Beneheart R1.<br/>ECG Trunk Cable and leads Set<br/>Mindray Beneheart R1</p>  |   |   |   |
| <p><b>ML7977</b><br/>Cable Holter <b>Spiderview</b> RC019<br/>ECG Holter Cable 3 LEADS<br/>XYZ / 7E (<b>7 Leads</b>) Ref. RC019</p>   |   | <p><b>ML7975</b><br/>RC019 cable for Holter<br/><b>Spiderview</b> Sorin<br/>ECG Holter Cable 3 LEADS<br/>XYZ / 7E (<b>5 Leads</b>) RC019</p>   |  |
| <p><b>RC93606</b><br/>Holter Cable 70869, 7-LEAD<br/>HOLTER MAIN CABLES, 0.60<br/>m. 70870, 1.20 m</p>  |   | <p><b>RC27606</b><br/>Din style holter 7-lead<br/>leadwires, <b>L=60 cm</b>, DIN 1.5<br/><b>RC27610</b><br/>Din style holter 7-lead<br/>leadwires, <b>L=100 cm</b>, DIN 1.5</p>  | <p><b>RC27912</b><br/>Din style<br/>holter 7-lead<br/>leadwires,<br/><b>L=120 cm</b>,<br/>DIN 1.5</p>  |





|  |  |   |   |
|--|--|---|---|
| <p><b>HT2793</b><br/>Holter 5-LEAD HOLTER MAIN CABLES, 0.60 m, 1.20 m<br/>Din style Holter 5-lead leadwires , L=60 cm , DIN 1.5<br/>Red, White, Green, Black and Yellow</p>  |   | <p><b>ML2156</b><br/>Holter 5-LEAD HOLTER MAIN CABLES, 0.60m 25452-513 (Holter)<br/>Red, White, Green, Black and Brown</p>  | <p>Din style<br/>Holter 5-lead leadwires , L=60 cm , DIN 1.5</p>                                      |
| <p><b>ML2289</b><br/>Holter, GE SEER MC Holter cable , 10-lead ,V1-V6, LL, LA, RA, RL, 80cm<br/><b>ML2389:</b> L=100cm <b>ML2394:</b> L=60cm</p>   |   |    | <p><b>ML2390</b><br/>Din style holter 10-lead leadwires , L=120cm , DIN1.5</p>                        |
| <p><b>ML6505</b><br/>Holter, Cable Part No. 17330 Del Mar Avionics for holter. Del Mar Avionics 17330, 17630. Cable Del Mar Avionics Part No. 17330 Without extension cordon. Del Mar Avionics Part No. 17630 With extension cordon 1,2m or 0,60m. ( I.P. 70632 or I.P. 70630)</p> |  <p>Del Mar Avionics Part No. 17630 With extension cordon 1,2m or 0,60m. ( I.P. 70632 or I.P. 70630)</p> |  <p>HOLTER 5 Voies. I.P. 70632 (1,2m) male 2 mm Jeu de 5 Pression IEC</p>                          | <p>Del Mar Avionics (Electrode &amp; Extension Cord Kit) Part No. 17630 for holter.</p>               |
| <p><b>ML2288</b><br/>Holter, 419598-005, GE SEER MC Holter cable , 5-lead 2 channel , 80cm</p>   |   |    | <p><b>ML2216</b><br/>Holter 5-LEAD HOLTER MAIN CABLES</p>   |
| <p><b>ML2287</b><br/>Holter, 419598-007, GE SEER MC Holter cable, 7-lead 3 Channel , 80cm</p>  |    |   | <p><b>RC27610</b><br/>Holter Cable 7-LEAD HOLTER MAIN CABLES, leadwires 1m</p>                        |
| <p><b>ML4596</b><br/>Cabo Holter para aparelho Holter 24H DigiTrak XT, Holter 5 lead cable, Philips Model 860322. Philips DigiTrak XT Patient Cables- 5 Leads - Medium Cable is standard 36 or 48 inches long, (89,44 or 121,92cm)</p>   |   | <p>Holter cable to device 24H DigiTrak XT Holter 5 lead cable, brand Philips Model: 860322.</p>  | <p>Philips/Zymed 3 channel 24hr DigiTrak- Plus &amp; Quinton Multitrack Digital Holter</p>            |
| <p><b>ML4144</b><br/>Cabo para registador Holter Schiller -2 canais MT100/3. Holter Schiller 2 canais MT-100/3 Schiller 5 lead patient cable IEC MT-100/3 Holter Shiller R6 2.400111</p>   |   |    | <p>Cabo para registador Holter Schiller MT100/3. Schiller 5 lead patient cable IEC MT100/3 Holter</p> |
| <p><b>ML4145</b><br/>Schiller MT200<br/>Schiller MT-200<br/>6 lead holter cable</p>  |   | <p><b>ML4146</b><br/>Schiller MT200<br/>Schiller MT-200<br/>4 lead holter cable</p>   |                  |
| <p><b>ML4613</b><br/>Cable MORTARA Holter 12 leadwires model HT0002. Cabo ECG Holter Mortara code 9293-017-51. MORTARA Holter H12, X12, X12+</p>   |   | <p>929301751</p>   | <p>Holter RS2011100156 5. Cabo Mortara 9293-017-51. 9293-017-51 - Set 10 leads, Snap Mortara</p>      |
| <p><b>ML4442</b><br/>Holter, Medilog Holter 7 lead cable , Redel 14pin -&gt;DIN /7</p>   |   |    | <p><b>RC27606</b><br/>Din style holter 7-lead leadwires , L=60 cm , DIN 1.5</p>                       |




|  |   |  |   |
|--|---|--|---|
| <p><b>ML1007</b><br/>Holter 2-lead wires, 2.5mm and 3.5mm plug. cable Holter Instromedix King of Hearts Express AF</p>   |   | <p>Instromedix<br/>Patient Cable-2 Lead<br/>#420-0073-00<br/>•For use with:•King of Hearts Monitor<br/>•King of Hearts Express•King of Hearts Express Plus<br/>•King of Hearts Express AF</p>  | <p>•Cable length: 30 inches<br/>•AHA standard<br/>•Dual phono plug 90 degree connector<br/>•Black/White leads</p> |
| <p><b>RC12288</b><br/>Holter, GE/Rozinn Holter 7 lead cable , Round 7pin -&gt;DIN/7 , 2 Channel ,80cm</p>  |    |   | <p><b>RC27606</b><br/>Din style holter 7-lead leadwires , L=60 cm , DIN 1.5</p>                                   |
| <p><b>RC22285</b><br/>Holter, GE/Rozinn Holter 5 lead cable , Round 7pin -&gt; DIN /5 , 2 channel, 80cm</p>  |    |   | <p><b>ML2216</b><br/>leadwires set,5-Lead ,Snap ,IEC ,40in ,DIN L=60 cm</p>                                       |
| <p><b>ML7162</b><br/>GE Seer 1000 20 Pin 7 Lead Holter Monitor Ecg Cable IEC AHA.<br/>GE: (Adult) 210413-001<br/>GE brand; Model: Seer1000; Cable Ref. GE: (Adult) 2104813-001.<br/>SEER 1000 Cable 2067634-001</p>  |    |   | <p>GE<br/>Seer1000 20 Pin 7 Lead Holter Cable 2067634-001<br/>Medke<br/>P/N: G71100S</p>                          |
| <p><b>ML5118</b><br/>Cable Holter DMS 300/3A<br/>Cabo Holter, Ref. SC-15F/7JA 20200213 DMS 300-3A<br/>Cabo Holter SC-15F/7JA<br/>DMS300-3A-7LCABLE DMS 300-3<br/><b>Holter 7 lead ECG leadwires,</b></p>   |  | <p><b>ML1615</b><br/><b>Holter 5 lead ECG leadwires</b></p>    | <p><b>ML1615</b><br/>DMS 300-3A DMS 300-3 Holter<br/><b>5 lead ECG leadwires</b></p>                              |
| <p><b>ML5117</b><br/>DMS 300-3A HDMI Holter Patient Recorder Cable, IEC. Snap. Holter DMS 300A. <b>SC-19F/7N, 20202516</b><br/>Holter DMS 303A, Holter DMS303A Holter DMS300A. , Holter DMS 304A<br/>DMS 300-3 Holter <b>7 lead ECG leadwires</b><br/>DMS-20202516</p> |  | <p><b>ML5177</b><br/>DMS 300-3A HDMI Holter Patient Cable Used with Holter recorders, applicable for the following models of Holter recorder: DMS300-3, 300-3A, Holter DMS 304A, Holter DMS 300. <b>DMS 300-3 -20202519</b><br/><b>Holter 5 lead ECG leadwires</b></p> |                              |
| <p><b>RL11905</b><br/>Holter, 2008594-003, GE SEER Light Patient Cable /Leadwires <b>2 Cannel (5-Ld) AHA.</b> GE Healthcare SEER Light Cable/Leadwires - Holter 2 Channel and 3 Channel Cables and Leadwires For GE SEER Light Holter Recorder (2008594-002)</p>       |  | <p>Description:<br/>•Dedicated cable/leadwire set eliminates chance for leadwires to be pulled away from cable<br/>•Lightweight cable/leadwire set for patient comfort<br/>•Available in 2 Channel and 3 Channel styles</p>  | <p>•SEER Light Holter Patient Cables/Lead wires<br/>•3 Channel (7-Ld), AHA</p>                                    |
| <p><b>RL11907</b><br/>Holter, 2008594-004 GE SEER Light Patient Cable /Leadwires <b>3 Cannel (7-Ld) AHA.</b> GE Healthcare SEER Light Cable/Leadwires - Holter 2 Channel and 3 Channel Cables and Leadwires For GE SEER Light Holter Recorder (2008594-002)</p>        |  | <p>Description:<br/>•Dedicated cable/leadwire set eliminates chance for leadwires to be pulled away from cable<br/>•Lightweight cable/leadwire set for patient comfort<br/>•Available in 2 Channel and 3 Channel styles</p>  | <p>•SEER Light Holter Patient Cables/Lead wires<br/>•3 Channel (7-Ld), AHA</p>                                    |
| <p><b>ML4914</b><br/>ECG Cable Holter <b>10 leads</b> for Biomedical Instruments de 3 channels BI9800TL+</p>   |  | <p><b>ML4914</b><br/>Cable Holter BI9800, BI9000, BI9800TL+, Holter ECG Cable 10 Lead Snap. Edan Holter ECG</p>  |                              |
| <p><b>ML5615</b><br/>Mortara 9293-044-51<br/>Cables for Mortara T12<br/>12 Lead Telemetry Module</p>   |  | <p>12 Lead Telemetry Module Only Cable</p>   |                              |




|  |  |  |  |
|--|--|--|--|
| <p><b>RL11900</b><br/>         Holter, BU010-1643-00, 7-leads set<br/> <b>7-lead</b> Patient Cable for Burdick Vision 5L holter recorder 3 Canais. BU010164300. Burdick Holter Leads Vision 5L holter recorder 3 Canais Burdick Vision 5L Holter Record Burdick Vision Holter patient Cable - <b>7 leads</b> Patient Cable 010-1643-00</p>   |    | <p>Cabo Holter para Burdick Vision Cable 010-1643-00 (BU010-1643-00)</p>    |   |
| <p><b>RL11901</b><br/>         Holter BU010-1642-00 5-leads set<br/> <b>5-lead</b> Patient cable for Burdick Vision 5L holter recorder 3 Canais. Burdick Holter Leads BU010164200. Spacelabs Vision 5L Vision 5L Digital Holter Record Burdick Vision Holter patient Cable - <b>5 leads</b></p>  |    | <p>Cabo Holter para Burdick Vision Cable 010-1642-00 (BU010-1642-00)</p>    | <p>10-1642-00 5-leads set. Patient Cable 010-1642-00</p>   |
| <p><b>RK3307</b><br/>         Holter Patient Cable 7- Leads 175-1546-00, 3 Channel, Safety, 92513 Holter, Burdick DIN Safty ECG, Trunk Cable, 7- Leads<br/>         Cable Holter 92513 Holter Burdick holter 7-lead patient fit for Din style LeadWirescable OEM p/n: 175-1546-00. HC030-7A. Burdick Holter Din 7-Ld Patient Cable, AHA, PCB 7pin to Din 5J. Burdick 92513 Holter Cable, ECG Din 7-Lead Patient Cable.<br/> <b>Cables with 45 and 80cm</b><br/> <b>Cabo para Burdick Holter 92513</b><br/>         Holter Spacelabs model 6732</p> |   |   | <p><b>RC27936</b><br/>         holter 7-lead leadwires ,<br/> <b>L=0,60 cm</b> ,<br/>         DIN 1.5<br/> <b>RC27610</b><br/>         holter 7-lead leadwires ,<br/> <b>L=100 cm</b> ,<br/>         DIN 1.5<br/> <b>RC27912</b><br/> <b>L=120 cm</b> ,<br/>         DIN 1.5</p> |
| <p><b>ML2001A</b><br/>         Welch Allyn Cable for ECG CP50 Welch Allyn<br/>         Welch Allyn CP 50 Welch Allyn<br/>         Welch Allyn CP-50 Welch Allyn<br/>         ECG cable 10 leads banana 4.0</p>   |  | <p>Cabo ECG CP 50 Welch Allyn Cabo ECG CP-50 Welch Allyn Cabo ECG CP50 Welch Allyn<br/> <b>10 leads banana 4.0</b></p>   |  |
| <p><b>ML1111</b><br/>         Welch Allyn ECG Cable and leadwires, IEC, <b>Banana</b> connector. ECG cable from Welch Allyn model SE-PRO-EUR-600, 10-Lead Banana. ECG cable from <b>Welch Allyn model SE-PRO-600</b>. (RE-PC-AHA-BAN 10-Lead Patient Resting Cable, AHA, 4 mm Banana for use with SE-PRO-600) <b>Redel 14 pin</b></p>  |  | <p><b>ML1111</b><br/>         Welch Allyn ECG Cable Model SE-PRO-EUR-600<br/>         Welch Allyn CABLE SE-PRO-600. <b>Redel 14pin</b><br/>         Cable BSK032C <b>10-Lead IEC/AHA Banana ECG Cable</b></p>  | <p>for use with SE-PRO-600</p>    |
| <p><b>ML1112</b><br/>         Welch Allyn ECG Stress Cable and leadwires, IEC, Clip connector. ECG cable from Welch Allyn model SE-PRO-EUR-600, 10-Lead Clip. ECG Cable 10-Lead Stress Patient Cable, AHA, Pinch Clip, for use with SE-PRO-600. ECG cable from <b>Welch Allyn model SE-PRO-600</b>. 14pin</p>  |  | <p><b>ML1112</b><br/>         Welch Allyn Cardio Perfect <b>Stress Patient Cable</b> for PRO-600 ECG Recorder, 10-Lead Stress Patient Cable (IEC/AHA), <b>CLIP</b> for use with SE-PRO-600, <b>Redel 14pin</b><br/>         Stress Patient Cable BSK032G</p> | <p>RE-PC-AHA-BAN 10-Lead Patient Resting Cable, AHA, <b>14pin Clip</b> for use with SE-PRO-600</p>   |
| <p><b>ML1114</b><br/>         Welch Allyn ECG Stress cable and leadwires, IEC, SNAP connector., ECG Cable 10-Lead Stress Patient Cable, AHA, Pinch SNAP, for use with SE-PRO-600. ECG cable from <b>Welch Allyn model SE-PRO-600</b>. RE-PC-AHA-BAN 10-Lead Patient Resting Cable, IEC/AHA, SNAP for use with SE-PRO-600</p>   |  | <p><b>ML1114</b><br/>         Welch Allyn Cardio Perfect <b>Stress Patient Cable</b> for PRO-600 ECG Recorder 10-Lead Stress Patient Cable (IEC/AHA), <b>SNAP</b> for use with SE-PRO-600<br/> <b>Redel 14pin</b></p>  | <p>RE-PC-AHA-BAN 10-Lead Patient Resting Cable, AHA, <b>14pin SNAP</b> for use with SE-PRO-600</p>   |
| <p><b>ML5294</b><br/>         Mortara H3+, Holter Cable<br/>         Mortara 9293-036-51.<br/>         H3+ 5-wire, 3-lead patient cable - IEC 27" (69cm). 9293-036-53.<br/> <b>ECG cable 5 Leads snap for Mortara H3, 69cm.</b></p>  |  | <p>Mortara 9293-036-53 H3+ 5-wire, 3-lead patient cable - IEC 27" (69cm). ECG cable 5 Leads snap for Mortara H3, 69cm.<br/> <b>9293-036-53</b><br/> <b>(Old code 9293-036-51)</b></p>  | <p><b>ML5294B</b><br/>         HT026 Mortara cable H3 leadwires. Holter cable P2800092</p>   |




|   |   |  |  |
|---|---|--|--|
| <p><b>ML5049</b><br/>Holter Ecg cable Ref 2.400161 from Holter Medilog AR12 Plus<br/><b>ECG cable 5 Leads snap for Medilog AR12 Plus.</b><br/>Schiller Plus 5 lead holter ECG cable for AR12 Plus AR4 Plus FD5</p>                          |   | <p><b>ML5047</b><br/>Holter Ecg cable from Holter Medilog AR12 Plus<br/><b>ECG cable 7 Leads snap</b><br/>Schiller Plus 7 lead holter ECG cable for AR4 Plus FD5</p>           |                             |
| <p><b>ML2444</b><br/><b>10 leads ECG cable Holter</b><br/>SCHILLER MEDILOG FD12 PLUS as picture attached. (cable (01) 07613365002430. (10)019453.</p>   |   | <p><b>ML2447</b><br/><b>7 leads ECG cable Holter</b><br/>SCHILLER MEDILOG FD12 PLUS as picture attached. (cable (01)07613365002447, (10)019184</p>                             | <p>HOLTER SCHILLER MEDILOG FD12</p>  |
| <p><b>ML2448</b><br/><b>4 leads ECG cable Holter</b><br/>SCHILLER MEDILOG FD12 PLUS as picture attached. (cable 2.400173</p>  |    | <p>Schiller 4-lead patient cable push button, FD12plus - 2.400173</p>  | <p>HOLTER SCHILLER MEDILOG FD12</p>  |
| <p><b>ML5684</b><br/>Welch Allyn USB Cable PRO 60023, PRO 60024, PRO 60025<br/>USB-Interface for Welch Allyn CardioPerfect.<br/>Welch Allyn PC-Based Exercise ECG USB Cable</p>   |    | <p>The USB-Interface for Welch Allyn Cardio Perfect is available in 3 different cable lengths:<br/>- PRO-60023, 2m<br/>- PRO-60024, 3m<br/>- PRO-60025, 5m</p>                 | <p>- PRO 60023<br/>- PRO 60024<br/>- PRO 60025</p>   |
| <p><b>ML3314</b><br/>FX-101, FUKUDA Denshi. Autocardiner FX-3010, FX-4010, FX-2111, FX-2155.<br/>FUKUDA EKG cable with 12 leads AHA or IEC DB15M&gt; <b>Din 3.0</b>(Original no: FX-101), Banana 4.0, Clip, Snap,<br/>See <b>ML3304</b></p> |   |   | <p><b>Din 3.0</b><br/>FCP-7411, FX-7402, FX-7202, FX-2111<br/>FX-7102, 603D / 603ER / 603S/ 606. I.P.71509</p> |
| <p><b>ML3704</b><br/>Fukuda Denshi<br/>FX-102F, FX-328UY, FCP 2155, FCP 2201, FCP 4101, FCP 4101A, CardiMax FX-7101, FX-7202, FX-7402. FX-702, FX 702<br/>Fukuda EKG cable , <b>AHA</b>, DB15M&gt; <b>Banana 4.0</b></p>                    |  | <p>Auto Cardiner Fukuda Denshi<br/>FX-102F, FX-328UY, FCP 2155, FCP 2201, FCP 4101, FCP 4101A, CardiMax FX-7101, FX-7202, FX-7402. FX-702, FX 702<br/>Fukuda CP-104T Cable</p> | <p>Fukuda EKG cable , <b>AHA</b>, DB15M&gt; <b>Banana 4.0</b></p>  |
| <p><b>ML1152</b><br/>One-piece Shield ECG cable, 3-Lead, <b>Snap, IEC</b>, 12ft, round <b>10pin</b>.<br/>ECG Cable for Drager ventilator monitor PM 8060 Vitara. PM8010 , PM8014, Babyguard 8000 , Cato, Babylog</p>                        |  |   | <p>8000/8000 plus PB8800 Series 8000, Cicero , Parameterbox UM1/3/3.1 , RM1</p>                                |
| <p><b>MC1284</b><br/>10-Lead ECG leadwires for beijing JincoMED holter, 11pin, snap, IEC</p>  |  |   |  |
| <p><b>ML4122</b><br/>Fukuda Denshi 5-lead ECG cable , IEC , Round 12pin<br/>Fukuda DYNASCOPE 7100, 7200, 7300. Beside Monitor DS-7101LT, DS-7141, DS-7101L</p>  |  |   | <p><b>ML4124</b><br/>Fukuda ECG leadwire set, 5-Lead, Grabber, IEC, 40IN &amp; 2pin</p>                        |
| <p><b>ML4129</b><br/>Fukuda Denshi Fukuda 3-lead trunk cable ,With simens 3-lead yoke , round 12pin -2PIN</p>   |  | <p><b>ML4130</b><br/>Fukuda 3-lead leadwires , snap , IEC type .Siemens leadwire connection</p>  |  |



|   |   |  |  |
|---|---|--|--|
| <p><b>ML2599</b><br/>Petas KM800 4-lead ECG cable with leadwires, IEC, R, N, F, L .Grabber, Redel 6pin with 1bit.<br/><b>One piece ECG cable</b></p>  |    |   | <p>PetasKM800<br/>Petas KM 800</p>   |
| <p><b>ML4018</b><br/>Mobimed EKG cable, IEC Round 12pin &gt; Banana 4.0. Defibrillator/Monitor</p>  |    |   | <p>Mobimed<br/>Defibrillator /Monitor</p>  |
| <p><b>ML4799</b><br/>Cable Mortara Quinton Q-Stress Mortara Q Stress EKG Cable 10 Lead Snap Cardiac Science Quinton EKG Cable 10 Lead Snap Mortara Quinton Q-Stress</p>   |    |   | <p>Mortara<br/>Quinton<br/>Q-Stress<br/>12 Pin</p>   |
| <p><b>RC6643</b><br/>Cabo ECG Nikon Kohden Cardiofax V Mod. 8240K (ECAPS 12). 72269, <b>MONOBLOC 10 LEADS BANANA, NIHON KOHDEN,</b> Cable ECG Nikon Kohden Cardiofax V ECAPS12, Mod. 8240K</p>  |    |   |                           |
| <p><b>ML4589</b><br/>Burdick EK-10<br/>Burdick EKG cable , AHA , DB15M&gt;Banana 4.0</p>  |   |    |  |
| <p><b>ML7497</b><br/>CONTEC ECG EKG Electrocardiograph for ECG Contec 300G. Contec Cable</p>   |   | <p>Cable 12-lead ECG. Cable CONTEC ECG EKG Electrocardiograph for ECG100G/ECG300G/ECG600G. Guardian ECG 300G, ECG-100G/ECG-300G/ECG-600G. ECG 100G/ECG 300G/ECG 600G. Contec K10R-NK2-BI . Caremed K10R-NK2-BI</p>               | <p>Contec<br/>300G</p>  |
| <p><b>ML4119</b><br/>Mortara Surveyor 10-lead EKG cable, IEC , Round 12pin -&gt;Snap<br/><b>ML4121</b><br/>Mortara Surveyor 10-lead EKG cableA ,HA ,Round 12pin &gt;banana 4.0</p>  |   |   | <p><b>ML4120</b><br/>Mortara Surveyor 10-Lead EKG cable ,IEC, Round 12pin -&gt;Grabber</p>                   |
| <p><b>ML4008</b><br/>SIEMENS 3368425<br/>SIEMENS 33 68 425<br/>ECG MONITORING PATIENT MEDICAL MONITOR CABLE 3368425<br/>SIEMENS DRAGER 33 68 425<br/>interface cable AS3368425</p>  |    |   | <p>Interface cable 3368425 for Monitor Siemens AS3368425 cable for SC 6000/9000 Siemens</p>                  |
| <p><b>ML3019</b><br/>EP00500 Equipotential bonding cable medical grade <b>length 5 m.</b> Highly flexible connecting cable with green yellow insulation for potential equalization according to IEC EN 60601-1, IEC EN 60364 7-710, IEC EN 11197, DIN 42801, DIN 42801-2.<br/><b>Length 500cm, Diameter 6 mm2</b></p> |    | <p><b>ML3018</b><br/>EP00250 Equipotential bonding cable medical grade <b>length 2,5 m.</b> Highly flexible connecting cable with green yellow insulation for potential equalization<br/><b>Length 250cm, Diameter 6 mm2</b></p> |                         |










|  |   |  |  |
|--|---|--|--|
| <p><b>ML3791</b><br/>Physio Control Lifepak 12/15<br/>ECG Cable from Defibrillator Physio<br/>Control Lifepak 12. Lifepak 12/15<br/><b>Complete ECG Cable</b><br/>(Including the trunk 4 lead ECG<br/>leadwires + 6 lead ECG)</p>  |    |   |   |
| <p><b>ML3711</b><br/>Physio Control Lifepak 12/15/20<br/>ECG Cable from Defibrillator<br/>Physio Control Lifepak 12.<br/><b>Complete ECG Cable</b><br/>(Including the trunk + 6 lead<br/>ECG leadwires + 4 lead ECG<br/>leadwire) 11111-000020<br/><b>IEC Leads are included</b></p> |    | <p>Physio Control Compatible<br/>ECG Trunk Cable -<br/>11111-000020<br/>Physio Control Lifepak 11,<br/>Lifepak 12, Lifepak 15 and<br/>Lifepak 20. For use with Trunk<br/>Cable TD-210270. 11110-<br/>000102 (4 Leads), 11110-<br/>000103 (6 Leads)</p> |   |
| <p><b>ML3710</b><br/>ECG Cable for Physio Control<br/>Lifepak 12 Defibrillator.<br/>Lifepak 12/15 Patient Cable for<br/><b>use with 4 or 12 lead</b><br/>attachments. MD-110270<br/>11110-000111, 11110-000110<br/><b>Leads are not included</b></p>                                 |    | <p>To Monitor:<br/><br/>To 4/6 Leads:<br/></p>                                   | <p>Lifepak12,<br/>Lifepak 12.<br/>Lifepak 12</p>   |
| <p><b>ML3714</b><br/>Physio Control Lifepak 12<br/><b>6 Wire Précordial Lead IEC</b><br/>6 lead ecg leadwire.<br/><b>6 lead 11110-000103.</b><br/><b>4 lead 11110-000102</b></p>   |   | <p><b>ML3712</b></p>   | <p><b>ML3712</b><br/>Physio Control<br/>Lifepak 12,<br/><b>4 Wire</b><br/><b>Précordial</b><br/><b>Lead IEC</b><br/><b>MDA4-120S</b></p> |
| <p><b>ML3715</b><br/>Physio Control Lifepak 12/15/20.<br/>ECG Cable 10 leads Lifepak 12.<br/>ECG/EKG Physio Control<br/>Lifepak12/15/120 10-lead<br/>12pins Cable Fit for Medtronic<br/><b>Complete ECG Cable</b></p>  |  |   |   |
| <p><b>ML0094</b><br/>ECG Cable RIB Omicron FT<br/>Surveyor<br/>5-lead ECG cable with leadwires ,<br/>Round 12pin , IEC, <b>Snap</b></p>  |  | <p><b>ML0094G</b><br/>ECG Cable RIB Omicron FT<br/>Surveyor<br/>5-lead ECG cable with<br/>leadwires , Round 12pin , IEC,<br/><b>Grabber</b></p>  |   |
| <p><b>ML1319</b><br/>Aquarius plus Adult finger clip<br/>sensor, 8 pin blue, L=3M.<br/>Northern Meditec Aquarius Plus<br/>Adult finger clip spo2 sensor<br/>8pin blue L=3M</p>   |  |   |   |
| <p><b>ML6204</b><br/>Grabber Electrode Adapter to<br/>connect ECG cable 4.0 banana<br/>plug. ECG converter kit 10<br/>Leads Banana&gt; Graber</p>  |  |   |   |
| <p><b>ML6203,</b><br/>Electrode Adapter, 4.0 Banana -<br/>&gt; Snap, L=10CM.<br/>ECG converter kit 10 Leads<br/>Banana&gt; Snap</p>  |  | <p><b>ML6202</b><br/>Electrode Adapter, 4.0<br/>Banana-&gt; Grabber,<br/>L=10CM.<br/>ECG converter kit 10 Leads<br/>Banana&gt; Graber</p>  |   |




|  |   |  |  |
|--|---|--|--|
| <p><b>ML7548</b><br/>Charmcare Prizm3 ECG cable and leadwire, IEC, with snap connector. Redel-4 Pin ECG Cable (3-lead)</p>   |    |   | <p>Charmcare Prizm 3 ECG Cable BET3356A</p>  |
| <p><b>RL12968</b><br/>HP one piece ECG cable with leadwires, IEC, 3-Lead, Grabber, 12pin.<br/>Philips 9898031 43171.<br/>One piece 3-lead cable, grabber</p>   |   | <p>Philips Medical System:<br/>78352C, 78354C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M 4735A.</p>  | <p>Philips 9898031431 71. One piece 3-lead cable, grabber</p>  |
| <p><b>ML6607</b><br/>HP one piece ECG cable with leadwires, IEC, 5-Lead, Grabber 12pin M8105A MP5 Monitor, M4735A heartStart XL defibrillator monitor, M3535A/M3536A HeartStart MRx Defibrillator Monitor</p>  |    | <p>Philips Medical System :<br/>78352C, 78354C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M 4735A, 863066, 863068, SureSigns VM6, VM8 ECG Modules, M1001A/B, M1002A/B,</p>                             | <p>M1002A/B Multi-measurement modules, M3000A, M3001A, M3002A MMS X2, M8102A MP2 Monitor</p>                   |
| <p><b>ML6608</b><br/>Philips Medical, one piece ECG cable with leadwires, IEC, 5-Lead, Snap, 12pin.<br/>M3000A, M3001A, M3002A MMS X2, M8102A MP2 Monitor, M8105A MP5. monitor Philips VM6</p>   |  | <p>Philips 78352C, 78354C, 78834C Legacy HP patient Monitors 862472 to 862475 C3 Monitor, 862478, 862479 C3 Monitor, 863064, 863065, 863066, 863068 SureSigns VM6, VM8 ECG Modules, M1001A/B, M1002A/B Multi-measurement modules</p> | <p>Monitor M4735A heartStart XL defibrillator monitor, M3535A/M3536 A HeartStart MRx Defibrillator Monitor</p> |
| <p><b>RL1117</b><br/>Philips MP70, HP M1669A 3-lead Trunk cable, AAMI &amp; IEC label included. Round 12PIN</p>  |  |   | <p><b>RL1119</b><br/>HP M1974A 3-lead leadwires, IEC, Snap.</p>  |
| <p><b>RL1115</b><br/>HP M1668A, HP-2585 HP 2585, 5-lead Trunk cable, AAMI &amp; IEC label included. Round 12PIN</p>  |  |   | <p><b>RL1114</b><br/>HP M1664A 5-lead leadwires, AHA, Snap</p>   |
| <p><b>RL1118</b><br/>M1972A, HP M1972A 3-lead leadwires, IEC, Grabber HP3-90P</p>  |  |   | <p><b>RL1114</b><br/>MC065B-5, HP M1664A 5-lead leadwires, AHA, Snap</p>                                       |
| <p><b>LM4499</b><br/>Schiller EKG Cable, IEC, DB15M&gt;Grabber. Suitable for Schiller Cardiovit: AT-102, CS-200, Cardiovit AT-10 plus, AT-104 PC, AT-110, Futuremed: All Cardiette: All Astron: All BTL: All, Del Mar: 910, 920, 930,940,942,950, 960,970, 980. Camina: All Cabo Prova de esforço p/ECG Cardiovit / Shiller CS-200, Shiller CS200 - Ergo Sipro</p> |  | <p>CABLE MONOBLOC 10 VOIES PINCE. SCHILLER 15 CONTACTS MALES IEC PROTEGE 10K + ECL.VIS COURTE. SCHILLER (ELECTROCARDIOGRAPHIE)</p>              | <p>Cabo ECG (72884). Sheller/ CARDIOVIT (AT-3 / AT-5 / AT-6 / AT-10 / AT-60 / CS-100 / CS-200)</p>             |





















|  |  |   |  |
|--|--|---|--|
| <p><b>LM4466</b><br/>Schiller EKG Cable, IEC, DB15M&gt;Snap. Suitable for Schiller Cardiovit : AT-102, CS-200, Cardiovit AT-10 plus , AT-104 PC, AT-110, Futuemed : All Cardiette : All, Astron : All, BTL : All<br/>Del Mar: 910, 920, 930,940,942,950, 960,970, 980. Camina: All<br/>Cable E10R-SH1-B/I., Shiller CS200</p>                        |   |    | <p><b>LM4499</b><br/>Shiller CS 200 - ERGO SIPRO Cardiovit / Shiller CS-200, Shiller CS200 - Ergo Sipro Shiller CS 200</p>   |
| <p><b>LM5548</b><br/>Schiller EKG cable, IEC/AHA, DB15M&gt;Banana 4.0. Suitable for Schiller Cardiovit: AT-102, CS-200, Cardiovit AT-10 plus, AT-104 PC, AT-110, Futuemed: All, Cardiette: All, Astron: All BTL: All, Del Mar: 910, 920, 930,940,942,950, 960,970, 980. Camina : All</p>   |   |    | <p>Cardiovit / Shiller CS-200, Shiller CS200 - Ergo Sipro Ergo Sipro E10R-SH1-B/I</p>  |
| <p><b>LM5544</b><br/>Schiller EKG cable, IEC, DB15M&gt;Banana 4.0, Schiller Cardiovit: AT-101, AT101.<br/>Schiller Cardiovit ref. F6725R AT1 Schiller AT-1, Schiller AT-2 Plus E10R - SH1- B/I, ECG Shiller Cardiovit AT-1 ref. 2400070</p>  |   | <p>Schiller AT-1, Schiller AT-2 Schiller AT-101 Schiller AT101 Schiller AT 101</p>  | <p><b>ML7706</b><br/><b>Pediatric suction Electrodes</b><br/>Metal diameter : 16mm</p>  |
| <p><b>ML7108</b><br/>Multi-use Limb Electrode for all types of ECG electrode, <b>Metal plating AgCl (Dark), IEC.</b><br/>For Multi-use, Suitable for all 3.0DIN, 4.0 Banana, Snap and Grabber, Limb Electrode<br/><b>(4pcs / set)</b></p>  |  <p><b>ML7190</b><br/>Pediatric Multi-use Limb Electrode with Screw, Plastic clamp Between upper and lower <b>(4pcs / set)</b></p> |  <p>For Multi-use, Suitable for all 3.0DIN, 4.0 Banana, Snap and Grabber<br/>Metal diameter : 23mm<br/>Ball diameter : 32mm</p>                            | <p><b>ML7707</b><br/><b>Adult Multi use</b> electrode , 3.0, 4.0 ,Snap , grabber ,Metal plating AgCl (Dark)<br/><b>6pcs / set</b></p>                                      |
| <p><b>ML3103</b><br/>Mortara Instrument: ELE 100, ELI 200, ELI 150 ELI100, ELI200, ELI150<br/>Mortara EKG cable, IEC, DB15M&gt; <b>Banana 4.0</b></p>  |   | <p><b>ML3132</b><br/>Mortara Instrument: ELE 100, ELI 200, ELI 150 ELI100, ELI200, ELI150<br/><b>MORTARA STRESS CABLE</b><br/>For X-Scribe II<br/>Mortara EKG CABLE and leadwires , IEC, <b>Snap</b></p>                                      |   |
| <p><b>ML3131</b><br/>MORTARA STRESS CABLE 10-lead with 4mm grabber leads ends. For X-Scribe II, X-Scribe DOS, Portrait<br/>Mortara EKG CABLE and leadwires , IEC, <b>Grabber</b></p>   |   |    | <p>Mortara EKG CABLE and leadwires , IEC, <b>Grabber</b></p>   |
| <p><b>ML3033</b><br/>Mortara ELI 250, ELI 350, 10 10-lead ekg leadwire (Para ECG Mortara ELI 150) part number 72228 (Mortara ref. 9293-032-51). (MORTARA - SERIE 150-250-400/ PROTEGE 10K-15 CONTACTS MALES<br/>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI 150, 250, 350, 10. ELI 10, ELI 250, ELI 350. EF142R</p>               |   | <p>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI150, ELI250, ELI350, ELI10. ELI 10, ELI 150, ELI 250, ELI 350. EF142R</p>  | <p>Mortara ELI250 10-lead EKG cable, IEC , DB15M - <b>&gt;Banana</b></p>   |
| <p><b>ML4622</b><br/>Mortara ELI 250, ELI 350, 10 10-lead ekg leadwire (Para ECG Mortara ELI 150) (MORTARA - SERIE 150-250-400/ PROTEGE 10K-15 CONTACTS MALES.<br/>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI 150, 250, 350, 10. ELI 10, ELI 250, ELI 350. <b>Stress cable</b><br/><b>Cabo prova esforço e/ou pediátrico</b></p> |   | <p>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI150, ELI250, ELI350, ELI10. ELI 10, ELI 150, ELI 250, ELI 350. EF142R<br/><b>Cabo prova esforço e/ou pediátrico</b><br/><b>Cable Stress test</b><br/><b>Stress cable</b></p> | <p>Mortara ELI250 10-lead EKG cable, IEC , DB15M - <b>&gt;Snap</b></p>   |




















|   |   |  |  |
|---|---|--|--|
| <p><b>ML4623</b><br/>Mortara ELI 250, ELI 350, 10 10-lead ekg leadwire (Para ECG Mortara ELI 150) (MORTARA - SERIE 150-250-400/ PROTEGE 10K-15 CONTACTS MALES<br/>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI 150, 250, 350, 10. ELI 10, ELI 250, ELI 350. <b>Stress cable</b><br/><b>Cabo prova esforço e/ou pediátrico</b></p>   |    | <p>Cabo ECG, 10 derivações para ECG Mortara Mod. ELI 150, for ELI150, ELI250, ELI350, ELI10. ELI 10, ELI 150, ELI 250, ELI 350. EF142R</p>  | <p>Mortara ELI250 10-lead EKG cable, IEC, DB15M - <b>&gt;Grabber</b><br/><b>Cable Stress test</b><br/><b>Stress cable</b></p>  |
| <p><b>LM4412</b><br/>Schiller EKG Cable, IEC, DB15M&gt;Grabber. Suitable for Schiller Cardiovit: AT-102, CS-200, Cardiovit AT-10 plus, AT-104 PC, AT-110, Futuremed: All Cardiette: All Astron: All BTL: All, Del Mar: 910, 920, 930,940,942,950, 960,970, 980. Camina: All</p>   |    |   | <p>Cabo ECG (72884). Schiller/ CARDIOVIT (AT-3 / AT-5 / AT-6 / AT-10 / AT-60 / CS-100 / CS-200)</p>  |
| <p><b>AM7733</b><br/>COLIN, BP 88 - 306 - 508<br/>Cabo ECG completo p/Monitor Colin Ref. 70612NGJ3 + Ref. 72310 de 3 derivações ou Ref. 70612NGJ5+ Ref. 72311 de 5 derivações), COLIN, BP 88 - 306 - 508, 6 MALE PINS, 72310 (3 Derivações) Colin ECG BP-88 Serie 300061 AY1004 N.7 ou 72311(5 Derivações) (AY0004) MAIN CABLE 3 LEADS WITHOUT LEADWIRES and MAIN CABLE 5 LEADS WITHOUT LEADWIRES</p>   |    |   | <p>72311 de 5 derivações<br/><b>Cabo</b><br/>p/Monitor Colin ECG BP-88 Série 300061 AY1004</p>   |
| <p><b>ML7409</b><br/>Colin Bp-306 / Bp88 3-lead Trunk Cable Use With Colin Bp88 / Bp306, Colin BP-306/BP88 3-lead trunk cable, connector: L-6pin/6pin, IEC, 9ft, Compatibility: Use with Colin BP88/BP306. Specifications: Trunk cable connector: type 6pin 1.5 plugs, Patient Leadwire material: 26AWG/1C+shield, Patient end terminal: 4.0 Snap, Compatible with COLIN BP-88, BP-306 BP-88S patient monitor. Monitor Collin BP88 ref. AG1002 nº 5 ( AG 1002 nº 5)</p> |  |   | <p><b>ML7411</b><br/>Colin ECG Leadwire set-Group, 5-LD, Grabber, IEC, 40in &amp; 6PIN<br/><b>ML7410</b><br/>Colin ECG Leadwire set-Group, 3-LD, Snap, IEC, 40in &amp; 6PIN</p>  |
| <p><b>ML7408</b><br/>Colin BP88S ECG Trunk cable, 5-Lead, IEC, 9ft, &amp;6pin plug/6pin yoke. (72330, AY1005 (TRUNK CABLE 3 LEADS WITHOUT LEADWIRES) and 72331, AY1005 (TRUNK CABLE 5 LEADS WITHOUT LEADWIRES): BP 88 NEW GENERATION 6 FEMALE PINS. Colin BP88, BP306<br/><b>Cabo Intermediário ECG 3 derivações para monitor Collin BP-S510</b></p>  |  |   | <p><b>ML7421</b>, Colin BP88, BP306, one piece ECG cable with <b>leadwire 3 -lead</b><br/><b>ML7413</b>, Colin BP88, BP306, one piece ECG cable with <b>leadwire 5 -lead</b></p> |
| <p><b>ML9456</b><br/>Colin BP88S, one piece ECG cable with leadwire 3-,lead,<b>IEC, Grabber</b>, Round 6in. <b>Cabo de ECG 3 derivações para monitor Collin BP-S510</b><br/><b>Direct connect</b></p>   |  | <p><b>ML9457</b><br/>Colin BP88S, one piece ECG cable with leadwire3-,lead, <b>IEC, SNAP</b>, Round 6i. <b>Cabo de ECG 3 derivações para monitor Collin BP-S510</b><br/><b>Direct connect</b></p>                              |   |
| <p><b>ML9458</b><br/>Colin BP88S, one piece ECG cable with leadwire, <b>5-lead,IEC, Grabber</b>, Round 6in. <b>Cabo de ECG 5 derivações para monitor Collin BP-S510</b><br/><b>Direct Connect</b></p>   |  | <p><b>ML9459</b><br/>Colin BP88S, one piece ECG cable with leadwire <b>5-,lead,IEC, Snap</b>, Round 6in. <b>Cabo ECG 5 derivações para monitor Collin BP-S510</b><br/><b>Direct Connect</b></p>                                |   |























|   |  |   |  |
|---|--|---|--|
| <p><b>RL12968</b><br/>Cabo para Monitor Philips SureSigns VM4, VM6, VM8 ECG Cables and Leads.<br/>ECG Philips 989803143171<br/>Cables monobloc 3 leads</p>  |     | <p>Philips/Agilent/HP, 12 male pins Protected, Cath Lab 9000 M 1275 A, M1275A Merlin 78300C, Viridia M3, Intellivue C3. 989803143171</p>  |  |
| <p><b>ML6607</b><br/>Cabo para Monitor Philips SureSigns VM4, VM6, VM8 ECG Cables and Leads.<br/>Cables monobloc 5 leads for monitor Philips VM6 ou VM4</p>   |     | <p>Philips/Agilent/HP, 12 male pins Protected, Cath Lab 9000 M 1275 A, M1275A Merlin 78300C, Viridia M3, Intellivue C3.</p>   | <p>CATHLAB 9000, M12 75A, Gamme Merlin Serie 78300 C 12 Male Contacts Protected</p>  |
| <p><b>ML8871</b><br/>Burdick Atria 6100<br/>Burdick Atria 6100 EKG cable , AHA, DB15M&gt;Banana 4.0<br/>Cabo para Electrocardiografo Burdick Atria 6100</p>   |     | <p>Cable for Electrocardiograph Burdick Atria 6100</p>  |  |
| <p><b>ML4824</b><br/><b>STORZ 495 NE</b> Fiberoptic Light Cable. Rigid Endoscope Light Source Fiber Optic Cable, 3m length, active diameter 4.8mm<br/><b>ML4819</b><br/><b>STORZ 495 NA</b> Fiberoptic Light Cable. The Storz cable 495 NA has an active diameter of 3.5mm and a length of 230 cm, Storz cable 495NA<br/><b>ML4820</b><br/><b>STORZ 495 NB</b> Fiberoptic Light cable. The Storz cable 495 NB has an active diameter of 4.8mm and a length of 180 cm. Storz cable 495NB</p> |    | <p><b>ML4588</b><br/>Storz 495 BD, Storz 495 BDM 12-297-ST<br/>Fiberoptic Light Cable</p>   | <p><b>ML4588</b><br/>Storz 495BD<br/>Storz 495BDM<br/>12-297-ST<br/>Fiberoptic Light Cable</p>  |
| <p><b>ML8807</b><br/>STORZ HOPKINS 8706CJ Endoscope<br/>STORZ Hopkins 8706CA 70 Strobo Rhinolaryngoscope<br/>8706 CJ Strobo-Laryngoscope with Hopkins Lateral Telescope 70° Working Length 17cm Oval Sheath, 7.2 mm x 9.3mm<br/>with intergral fiber optic light cable</p>  |   | <br>  | <p>KARL STORZ HOPKINS 8706CJ Endoscope</p>   |
| <p><b>ML0202</b><br/>Cold light cable 104-252-000, Item number: 104-252-000.<br/>For Kaps SOM 42/52/62. Microscopes. Incl. kink protection.<br/><b>Diam. = 6,80mm. (L = 1400 mm)</b><br/>Kaps 104-252-000</p>   |   | <br>  | <p>Cold light cable<br/><b>(L = 1400 mm)</b><br/><b>Kaps 104-252-000</b></p>   |
| <p><b>ML1108</b><br/>Cabo Bisturi eléctrico Berchtold Eletrobom-Integra EZ 294-45. EZ24945<br/>Cabo Bisturi p/ligar Placa descartável ref. T210 Commed. P/N 40421A. CABLE POUR PLAQUE BISTOURI. JACK 6.35mm<br/>BERCHTOLD BOVIE LAMIDEY MARTIN ERBE A.G. L3M. Reference Constructeur: GN249. Type: GK055 - GN060 - GN300. Marque: Aesculap. Commed EZ 294-45. Adapter cable - 3m, Cabo de placa para Bisturi Birtcher 5000.</p>   |   | <p><b>ML1102</b><br/>Disposable Electrosurgical-Monopolar Patient Plate,<br/>Size: <b>104*195mm (Adult)</b><br/><b>91*145mm (Pediatric)</b></p>  | <p>Use with Cable <b>ML1108</b></p>  |
| <p><b>ML1109</b><br/>Cabo Bisturi p/ligar Placa Neutra Commed 180 ref T229. (Cabos para bisturi eléctrico Valleylab T229 ) ELECTROBISTURI 98EGE088, EXCALIBUR PLUS</p>  |  | <p>40427A. Cabo Bisturi p/ligar Placa . CABLE POUR PLAQUE BISTOURI. JACK 2,5mm. Comprimento 3mt. Para Conmed180. EP01, P/N 40427A de 3 metros</p>   | <p>Tyco Ref. E 0560, P/N 40427A de 3 metros</p>  |




|  |   |   |  |
|--|---|---|--|
| <p>ERBE 20132-158 for APC applicators with APC connecting cable 20132-158</p>  |    | <p>Electrosurgery and Argon-Plasma Coagulation <b>ERBE</b></p>  |  |
| <p>ERBE 20194-077 (4m) e 20194-078 (5m).VIO, ICC, ACC, T-series, standard, for patient plates with contact stud. VIO, ICC, ACC, T-Serie, Standard, für Neutralelektroden mit Kontaktzunge 4 m 20194-077</p>  |    |    |  |
| <p>ERBE 20194-081 (4m), Martin standard, for patient plates with contact stud. Martin Standard, für Neutralelektroden mit Kontaktzunge 4 m 20194-081</p>   |    |    |  |
| <p>ERBE 20194-079 (4m) e 20194-086 (5m) para VIO, ICC 30 (ICC30), ACC, non-ERBE units, International, for patient plates with single contact surface and contact stud. VIO, ICC, ACC, Fremdgeräte, International, für ungeteilte Neutralelektroden mit Kontaktzunge 4m 20194-079</p>                                     |    |    |  |
| <p><b>ML9129</b> (4m) <b>ML9130</b> (5m)<br/>ERBE 20194-080 (4m) and ERBE 20194-087 (5m), VIO, ICC, ACC, non-ERBE units, international, for patient plates with split contact surface and contact stud. VIO, ICC, ACC, Fremdgeräte, International, für zweigeteilte Neutralelektroden mit Kontaktzunge 4 m 20194-080</p> | <br> | <br> |  |
| <p>ERBE 20194-068 (5m). VIO, ICC, ACC, standard, for patient plates with contact stud, with ECG connection . VIO, ICC, ACC, Standard, für Neutralelektroden mit Kontaktzunge, mit EKG-Anschluss 5 m 20194-068</p>  |   |   |  |
| <p>ERBE 20192-113. Monopolar connecting cable, VIO, ICC, ACC, international (Bovie Jack). for CUT and COAG, MIS instruments, length 4 m. max 5 kVp . 21960122</p>  |   |   | <p>85100-048<br/>max 95°C.<br/>30192-126<br/>max 138°C</p>                           |
| <p>350-040 Bipolar connecting cable (Compatible with Erbe ICC International &amp; US standard generators) 4.5 metre length. Logic cat. n: 610-2074, P/N: 350-040<br/>Description: cable for Erbe ICC 4,5mt<br/>Suitable on: Spinotto 4mm / Attacco europeo 80100011.</p>   |    |    |  |
| <p><b>RC6720</b><br/>CABLE ADAPTERS FOR MONOPOLAR HF SURGICAL INSTRUMENTS, European standards. Autoclavable (134° C / 18 mn). These cables are built with a 4mm female socket protected with a retractable sleeve, giving perfect security. No metallic part should</p>  |    | <p>Cabo monopolar para Laparoscopia ref.E0010R Ethicon. Cabo 1848.<br/>Cabo bisturi p/ toroscopia Ethicon E0010R. BOVIE VALLEYLAB</p>   |  |




|   |   |  |   |
|---|---|--|---|
| <p><b>ML4755</b><br/>Cabo Bisturi p/ligar Pinça<br/>Bipolar Bisturi Birtcher 5000<br/>(Commed). Ref. 37 01 35K (ref.<br/>370135K Bipolar Cable, 4,5m).<br/>HF 9450-50, HF9450-50.<br/>1441212F (5m). 1441211F (3m)</p>  |    |   | <p>Cabo Bisturi<br/>p/ligar Pinça<br/>Bipolar Bisturi<br/>Birtcher 5000<br/>(Commed BP-<br/>3500). For use<br/>with BP-3000.</p>              |
| <p><b>ML9351</b><br/>Bipolar connection cable, 5.00m<br/>for VALLEYLAB / ERBE<br/>80100013.<br/>Bipolar Cable Bissinger<br/>80100013</p>  |    |   |   |
| <p><b>ML2118</b><br/>CABLE ADAPTERS FOR<br/>MONOPOLAR HF SURGICAL<br/>INSTRUMENTS.<br/>BOVIE VALLEYLAB</p> <p>-MALE PLUG: Diameter: 8mm<br/>-FEMALE SOCKET Diameter: 4mm<br/>-CABLE LENGTH: 3,5m</p>  |    | <p>European standards.<br/>Autoclavable (134° C / 18 mn).<br/>These cables are built with a<br/>4mm female socket protected<br/>with a retractable sleeve,<br/>giving perfect security. No<br/>metallic part should be in<br/>contact with the user,<br/>according to current directives</p> | <p>Cabo<br/>monopolar<br/>para bisturi<br/>eléctrico<br/>(ficha macho<br/>8mm, ficha<br/>fêmea 4mm)<br/>Cabo com<br/>3,5m<br/>comprimento</p> |
| <p><b>ML7470</b><br/>Bipolar cable P283-284-10.<br/>Adapters for Forceps - European Standard<br/>Connection, LENGTH: 4.50 m, Fixed<br/>distance 29 mm, GENERATOR TYPE,<br/>EMC DOLLEY, LAMIDEY, VALLEYLAB,<br/>BOWA, 40493IP. Description: L=10 ft. (3<br/>m) 2"Banana Plug to Euro style con.<br/>Fit Euro style bipolar forceps</p> |   |    | <p><b>ML7470</b><br/>P283-284-10<br/>Bipolar<br/>Electrode<br/>Cable.</p>   |
| <p>STORZ N° 26176LA.<br/>No 700-351-245, Banana Pin<br/>male ø 4mm (Valleylab)<br/>(RZ700351245 ) Bipolar Cable<br/>4,5mm Valley. Pour pince<br/>bipolaire de Laparoscopie.<br/>STORZ-Kleppinger female 5<br/>mètres</p>  | <p>Pour pince bipolaire de Laparoscopie</p>  <p>STORZ-Kleppinger female 5 mètres</p> | <p>Banana Pin male ø 4mm<br/>(Valleylab)</p>    | <p>(STORZ N°<br/>26176LV)<br/>No 700-351-<br/>250</p>    |
| <p>No 700-350-050 , Europeans<br/>«Flat Plug» femelle 5 mètres</p>  |  <p>Europeen «Flat Plug» femelle 5 mètres</p>  | <p>Banana Pin male ø 4mm<br/>(Valleylab) No 700-350-050</p>   | <p>No 700-351-<br/>151(29mm)</p>                         |
| <p>No 700-351-045, Pour pince<br/>avec 2 connecteurs (PIN)</p>  | <p>Pour pince avec 2 connecteurs (PIN)</p>  <p>Droit 5 mètres</p>                    | <p>Banana Pin male ø 4mm<br/>(Valleylab) No 700-351-045</p>   | <p>No 700-351-<br/>151(29mm)</p>                         |
| <p>No 700-345-405, Pour pince<br/>avec 2 connecteurs (PIN)</p>  | <p>Pour pince avec 2 connecteurs (PIN)</p>  <p>Coudé 45° 5 mètres</p>                | <p>Banana Pin male ø 4mm<br/>(Valleylab) No 700-345-405</p>   | <p>No 700-351-<br/>251 (29mm)</p>                        |



|  |  |   |   |
|--|--|---|---|
| <p>Pour instruments de laparoscopie STORZ - Olympus - ACMI - etc. Connecteur type femelle ø 4mm 3 mètres</p>   | <p>Pour instruments de laparoscopie STORZ - Olympus - ACMI - etc.</p>  <p>Connecteur type femelle ø4mm 3 mètres</p> | <p>No 700-360-030, Adaptateur Bovie ø 8 mm (Valleylab) Adaptateur Bovie ø 8 mm (Valleylab) (STORZ N°26006M)</p>  <p>Adaptateur Bovie ø 8 mm (Valleylab) (STORZ N°26006M)</p> |   |
| <p>Cabo 10 derivações p/mesa de Angiografia Drager. Siemens ref. 6623818 E 295E ( 66 23 818 E295E) Nova ref. MP00780 (NORMA AMERICANA (AHA) (IEC2)</p>   |    | <p>Cabo 10 derivações p/mesa de Angiografia Drager. Siemens ref. 6623784 E 295E (66 23 784 E295E) Nova ref. MP00756 (NORMA EUROPEIA (IEC) (IEC1)</p>  |   |
| <p>Cabo ECG Précordial Siemens ref. 6627504 E295E, NOVA REF. MP00765 (NORMA EUROPEIA (IEC) (IEC1), para Cabo 10 derivações p/mesa de Angiografia Drager. Siemens ref. 6623818 E 295E (66 23 818 E295E) (NORMA EUROPEIA (IEC) (IEC1) OU (NORMA AMERICANA (AHA) (IEC2)</p> |   |    | <p>Interface Cable Siemens 6627504</p>          |
| <p>Cabo ECG 3 derivações DRAGER MS14555 MS 14555 SC 5000 / 6000 / 6000 P / 6002 / 7000 / 8000 / 9000 Delta / Delta XL, Gamma / Gamma XL / Gamma X XL, Kappa / Kappa XLT, Vista / Vista XL</p>  |   |   | <p>Interface Cable DRAGER MS14555</p>           |
| <p>2001292-002 E9007GN Multi-Link 3-Lead ECG Cable with Integrated Grabber Leadwires, IEC, Multi-Link Select, 12 ft. (3.6M)</p>  |   | <p>Cabo ECG refª E9007GN 2001292-002, para monitor marca Multiparametro GE, modelo Medical DASH 2000/2500</p>   |   |
| <p>Cabo Bisturi p/ligar Placa Neutra Commed 180. (Cabos de placa para bisturi eléctrico Valleylab)</p>   |   | <p>40427A. Cabo Bisturi p/ligar Placa . CABLE POUR PLAQUE BISTOURI. JACK 2,5mm. Comprimento 3mt. Para Conmed180. EP01, P/N 40427A de 3 metros</p>   | <p>Tyco Ref. E 0560, P/N 40427A de 3 metros</p> |
| <p><b>ML1902</b><br/>GE CO: Cardiac Output Cables, GE Medical, 2025248-002, Cardiac Output <b>Cable 2.7m</b></p>   |   | <p><b>ML1903</b><br/>GE CO: Cardiac Output Cables, GE Medical, 2025248-002, EC043B, Cardiac Output <b>Cable 3.6m</b></p>  |   |
| <p><b>ML1904</b><br/>Spacelab Medical, ONE PIECE EKG Cable With leadwires, AHA, Round 17pin&gt;<b>Snap</b></p>   |   | <p><b>ML1905</b><br/>Spacelab Medical, ONE PIECE EKG cable With leadwires, AHA , Round 17pin &gt;<b>Banana 4.0</b></p>  |   |




|   |   |  |  |
|---|---|--|--|
| <p><b>ML2243</b><br/>545300, 3-wire set of patient leads for above trunk cable, With clamps, 75 cm long, shielded, color coded: red, green, yellow, Original GE (Datex) # 545 300.<br/>Datex-Ohmeda AS/3, CS/3 Cardiocap, Saltite plus S/5 light monitor. <b>For IEC</b> ask <b>ML2243E</b><br/><b>IEC Color:</b> Red, Green and Yellow<br/><b>AHA Color:</b> Black, Red and White<br/>Cable 545300.<br/><b>RL4425</b><br/>3-wire set of patient leads, <b>Snap</b></p> |    |  <p><b>RL4425</b><br/>545315, 3-wire set of patient leads, <b>grabber</b> ref: 545315, monitores DATEX-ESGSTRON AS 3 e CARDIOK II.<br/><b>RL4426</b><br/>3-wire set of patient leads, <b>IEC Snap</b></p> | <p><b>RL4425</b><br/>545315, 3-wire set of patient leads, <b>grabber</b> ref: 545315, monitores DATEX-ESGSTRON AS 3 e CARDIOK II.<br/><b>RL4426</b><br/>3-wire set of patient leads, <b>Snap</b></p> |
| <p><b>ML4416</b><br/>One piece ECG cable with leadwire, <b>3-lead, IEC, Grabber</b>, 10 pin female. Original GE (Datex) Datex-Ohmeda AS/3, CS/3, Cardiocap, Saltite plus S/5 light monitor. <b>Direct connect</b></p>   |    | <p><b>IEC Color:</b><br/>Red, Green and Yellow<br/><b>AHA Color:</b><br/>Black, Red and White</p>  |  |
| <p><b>RH3342</b><br/><b>545306, 545301</b>, Datex-Ohmeda <b>5-Lead Trunk Cable, IEC</b>.<br/>Datex-Ohmeda AS/3, CS/3 Cardiocap. Saltite plus S/5 light monitor. Datex-Ohmeda 300-Series <b>545305</b>, Datex-Ohmeda 5-Lead Trunk Cable, AHA. Ref. <b>ML4256</b></p>                               |   |    | <p><b>RL4423</b><br/>545316, Datex Pro1000 leadwire set, <b>5-Lead, grabber, IEC</b>, 100CM and 2PIN. LW-36000MX/5I</p>  |
| <p><b>ML4417</b><br/>One piece ECG cable with leadwire, <b>5-lead, IEC, Grabber</b>, 10 pin female<br/>Datex-Ohmeda AS/3, CS/3, Cardiocap, Saltite plus S/5 light monitor. <b>Direct connect</b></p>  |  |   | <p><b>ML4410</b><br/>One piece ECG cable with leadwire, <b>5-lead, AHA, Snap</b> 10 pin female <b>Direct connect</b></p>   |
| <p><b>RH3304</b><br/>72002 (545301 or 545909) 5 LEADS MAIN CABLE WITHOUT LEADWIRE, DATEX, CARDIOCAP-AS 3-CS 3<br/>Datex-Ohmeda 5-Lead Trunk Cable, IEC.</p>   |  |   | <p><b>ML2391</b><br/>Cabo ECG Précordial de 5 cabos clip<br/>Référence Constructeur 5956466</p>  |
| <p><b>ML5313</b><br/>Cabo VGA Macho/Macho 15 pinos para Datashow Toshiba TDP-S80 , de 1,8m até 20metros, normalmente usa de 1,8m</p>  |  |  |  |
| <p>Cabo do Marcador de Eventos Cardiacall Delmar Reynolds ref. 23-0059 , 1 Channel, 2 lead patient loop cable CardioCall Individual Item</p>  |  |   |  |






|  |  |   |  |
|--|--|---|--|
| <p>Cabo Marcador de Eventos 23-0051- CABO 3 ELÉCTRODOS P/ Cardiocal, Cabo Marcador de Eventos, 2 canais 23-0051- Marcador eventos Spacelabs Healthcare.</p>            |  <p>SPACELABS<br/>Healthcare<br/>PATIENT CABLE 2 CHANNEL,<br/>PRODUCT CODE 23-0051<br/>Use only with CardioCall</p> |    |   |
| <p><b>ML9917</b><br/>Cable 5 Electrodes for Faros 360.<br/><b>Bittium Faros 360</b><br/>Ref. 900867<br/><b>ML9917A</b><br/>Cable 3 Electrodes for Faros Ref.900866</p> |    |    | <p><b>USB Cable 900826</b></p>   |
| <p><b>ML8717</b><br/>Cabo H.P 14307B c/8 pinos e 4 Leads clips +- 60cm terminal with 2mm Diameter. Cores: black, green, withe and grey. Cable 14307B</p>               |   |    | <p><b>ML7063</b><br/>Leadwire set (4) for Cable 14307B</p>   |
| <p><b>ML6602</b><br/>NEC 47524 3-lead Trunk Cable. NEC San-Ei Defb.model cardiopac 3m33. 47524 Trunk cable+ 47504 leadwire</p>   |    |   | <p><b>ML6603, NEC 47504</b><br/>3-lead IEC leadwires,</p>  |
| <p><b>ML6604</b><br/>NEC YCE205<br/>3-lead IEC Snap leadwires</p>  |   |   |  |
| <p><b>ML5504</b><br/>Datascope -Mindray, Passport V 0012-00-1745-01<br/>5-lead Trunk cable , Rectangular 12pin , 0012-00-1745-01,</p>                                  |   |  |  |
| <p><b>ML5515</b><br/>Datascope -Mindray, Passport V 0012-00-1503-02, 5-lead leadwires , AHA, Snap , use with <b>ML5504</b></p>   |   |  | <p><b>ML5516,</b><br/>Datascope - Mindray, Passport V 0012-00-1503-02, 5-lead leadwires , IEC, Snap , use with <b>ML5504</b></p> |
| <p><b>ML2915</b><br/>Mortara 10 lead leadwires, Banana 4.0, IEC, 10pcs/set ELI 230, Mortara ELI 280. Mortara 9293-046-61 Cable. Mortara WAM2500 Mortara WAM 2500</p>   |    |  | <p>ELI 150C, ELI 230, ELI 250C ELI 280, ELI 350</p>  |


























|   |   |  |   |
|---|---|--|---|
| <p><b>ML5517</b><br/>Datascope -Mindray, Passport V<br/>One-piece Shield ECG cable<br/>5-Lead ,Grabber ,IEC, 12ft<br/>Rectangular 12PIN</p>   |                                      |   | <p><b>ML5518</b><br/>Datascope -<br/>Mindray,<br/>Passport V<br/>One-piece<br/>Shield ECG<br/>cable, 5-Lead<br/>Snap ,IEC, 12ft<br/>12PIN</p> |
| <p><b>ML4814</b><br/>Mindray, Beneview T5. T6. T8<br/>One-piece Shield ECG cable,<br/>5-Lead, Grabber, IEC, 12ft,<br/>round 12PIN. Mindray,<br/>Beneview T5. T6. T8, PM5000,<br/>PM 6000.</p>   |                                      | <p>Cabo Monitorização<br/>Cardíaca<br/>12 Pins – 5 Derivações<br/>P/N 0010-30-42179.</p>   | <p>Mindray,<br/>Beneview<br/>T5. T6. T8,<br/>PM5000,<br/>PM 6000,</p>   |
| <p><b>ML4885</b><br/>Mindray, Beneview T5. T6. T8<br/>One-piece Shield ECG cable, 5-<br/>Lead, Snap, IEC, 12ft, round<br/>12PIN. Mindray, Beneview T5.<br/>T6. T8, PM5000, PM 6000,</p>   |                                      | <p>Cabo Monitorização<br/>Cardíaca<br/>12 Pins – 5 Derivações<br/>P/N 0010-30-42179.</p>   | <p>Mindray,<br/>Beneview<br/>T5. T6. T8,<br/>PM5000,<br/>PM 6000,</p>   |
| <p><b>ML7121</b><br/>Mindray 10lead EKG cable,<br/><b>IEC, Din 3.0.</b><br/>Mindray DECG-03A, DECG-030A.<br/>Spacelabs SL6, Edan SE-3, SE-12, SE-12<br/>Express, SE-601A, SE-601B, SE-601C, SE-<br/>1010PC, SE-1200, SE-12 Express, ST-1212<br/>and Nihon Kohden BSM-2301K.</p> |  <p><b>ML7121</b> <b>ML7122</b></p> | <p><b>ML7122</b><br/>Mindray 10lead EKG cable,<br/><b>IEC, Banana 4.0.</b><br/>Mindray DECG-03A, DECG-030A.<br/>Spacelabs SL6, Edan SE-3, SE-12,<br/>SE-12 Express, SE-601A, SE-601B,<br/>SE-601C, SE-1010PC, SE-1200, SE-<br/>12 Express, ST-1212 and Nihon<br/>Kohden BSM-2301K.</p> |   |
| <p><b>ML8880</b><br/>Mindray, Beneview T5. T6. T8<br/><b>5-lead Trunk cable</b>, 12pin Green<br/>connector, IEC/AHA included.<br/>Mindray, Beneview T5. T6. T8,<br/>PM5000, PM 6000.</p>  |                                    | <p><b>ML8883</b><br/>Mindray, Beneview T5. T6. T8<br/><b>3-lead Trunk cable</b>, 12pin<br/>Green connector, IEC/AHA<br/>included. Mindray,<br/>Beneview T5. T6. T8,<br/>PM5000, PM 6000.</p>   |   |
| <p><b>ML3633</b><br/>Mindray ONE PIECE EKG cable<br/>With leadwires, IEC, Round<br/>12pin &gt;Banana 4.0</p>  |                                    |   | <p><b>RL4424</b><br/>PRO 1000<br/>3-Lead ECG<br/>cable , IEC,<br/>12ft , &amp; 6PIN<br/>PRO1000</p>   |
| <p><b>ML4015</b><br/>SCHILLER ARGUS LCM<br/>Cabo Intermédio de um monitor<br/>SCHILLER ARGUS LCM.<br/>MAIN CABLE 3 LEADS 12 pins<br/>WITHOUT LEADWIRE</p>   |                                    |   | <p><b>LM8822</b><br/>Leadwires<br/>Clip</p>   |
| <p><b>ML5210</b><br/>SCHILLER ARGUS LCM<br/>Cabo Intermédio de um monitor<br/>SCHILLER ARGUS LCM.<br/>MAIN CABLE 10 LEADS WITHOUT<br/>LEADWIRE<br/>Schiller SCOPE CM8 – TM7, 12<br/>MALE PINS PROTECTED<br/>Schiller SCOPE CM8-TM7</p>  |                                    |   | <p><b>ML2488</b><br/>Leadwires<br/>Banana</p> <p><b>ML4419</b><br/>Leadwires<br/>Clip</p>   |




|  |  |   |   |
|--|--|---|---|
| <p><b>ML5542</b><br/>SCHILLER ARGUS LCM.<br/>3-lead ECG cable with leadwire ,<br/>12pin , IEC , Clip<br/><b>one piece Cable with Leadwires</b><br/>Schiller scope CMB-TM7</p>  |    |    | <p><b>ML5221</b><br/>SCHILLER ARGUS LCM.<br/>3-lead ECG cable with leadwire , 12pin , IEC , Snap<br/>Schiller scope CMB-TM7</p> |
| <p><b>ML1229</b><br/>H.P. M1600A. 5-lead ECG patient trunk cable for Philips monitors with 8-pin ECG input connector, 9' (2.7m) long.<br/>PHILIPS AGILENT H-P F/ SERIE43100/78300/78500/78600 DEF.<br/>CODMASTER PROTECTED</p>   |    | <p>Use with Philips Instrument:<br/>783XX series, 7835XA/B,<br/>7853X series, M1722A/B,<br/>M1723A/B M3500B, M5500B.<br/>ECG Trunk Cable<br/>M150A/M1600A.<br/>ECG Trunk Cable<br/>M1580A/M1600A</p>                                  | <p>For Philips instrument:<br/>M3500B (Heartstream/HeartStart XLT),<br/>M5500B (Heartstart 4000).</p>                           |
| <p><b>ML9995</b><br/>Philips ECG 5 lead Trunk cable M1530A<br/>Philips medical instruments:<br/>M1001A/B, M1002A/B, M1722A/B,<br/>M1723A/B, M2475B, M3000A,<br/>M3001A, M4735A, 78352C, 7834C,<br/>78834C<br/>See ML2002</p>   |     |    | <p>See ML1002</p>   |
| <p><b>ML5015</b><br/>Varidia, Merlin, 24C, V24CT, M1204A, 77010AC, 77020A/AC, 77020AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 78331A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A, 78353A/B, 78354A, 78352A/B, 78534A/B/C, M1722A/B, M1723A/B, M2400A, M2406A, M3500B &amp; More.</p>  |    | <p>12-PIN CABLE WITH LEADS<br/>~FITS MANY ~ VIRIDIA 24C<br/>M1204A Trunk Cable for ECG<br/>78346A, 78351T, 78352A,<br/>78353A/B, 78354A, 78352A/B,<br/>78534A/B/C, M1722A/B,<br/>M1723A/B, M2400A, M2406A,<br/>M3500B &amp; More.</p> |   |
| <p><b>ML5016</b><br/>Philips 8 pin Connector: Merlin, Viridia, Omnicare CMA, Heartstream, Codemaster XL, 77000AC, 77020A/AC, 77020AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 78331A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A, 78353A/B, 78354A, 78532A/B, 78534A/B/C, M1722A/B, M1723A/B, M2400A, M2406A, M2475B, M3500B</p> |  |    |    |
| <p><b>ML1012</b><br/>Cabo para HP 78353A Description:<br/>HP M1634A IBP Adapter cable,<br/>Hp12pin hp 4pin. Part No.:<br/>IBP001E. Cabo linha arterial (Ref. M1634A). Cabo Pressões p/monitor Philips Hewlet Packard (Cabo Ponta Vermelha).</p>  |   | <p>Part No.: IBP001E.<br/>Cabo linha arterial<br/>Ref. M1634A. Cabo Pressões p/monitor Philips Hewlet Packard (Cabo Ponta Vermelha).</p>  |   |
| <p><b>ML5537</b><br/>Spacelabs Medical, OEM Part No.700-0008-09. OEM Model: 90367, 90369, 90492, 90496, 91496, Ultraview. 3-Lead Tru-Link® Ultraview comp. Shielded ECG Cable (detachable leads) R, N, L.; Length 10ft. / 3.0m. ,PART NUMBER: CB-61396R/I</p>  |   | <p>Cable (detachable leads) R, N, L.; Length 10ft. / 3.0m. ,PART NUMBER: CB-61396R/I</p>  |   |




|  |  |   |  |
|--|--|---|--|
| <p><b>ML5600</b><br/>700-0008-07 Tru-Link ECG cable , 5-Lead , IEC/AHA included, 10ft &amp; 17PIN, MED, Spacelabs SL3-90S-I , 5 Lead wire snap. (Lead wire length: 0,9 m) Spacelabs: 90701, 90303, 90305, 90308, 90309, 90385, 90496 Ultraview, 90369 Ultraview, 90367 Ultraview</p> |   <p><b>ML5602</b></p> |    | <p><b>ML5601</b><br/>700-0007-09 Tru-Link Leadwire sets , 5-Lead , Snap , IEC, 24in &amp; 2PIN</p>       |
| <p><b>ML5611</b><br/>Tru-Link ECG cable , 3-Lead , IEC/AHA included, 10ft &amp; 17PIN, MED, Spacelabs SL3-90S-I , 3 Lead wire snap. (Lead wire length: 0,9 m) Spacelabs: 90701, 90303, 90305, 90308, 90309, 90385, 90496 , Ultraview, 90369 Ultraview, 90367 Ultraview</p>           |   <p><b>ML5722</b></p> |    | <p><b>ML5725</b><br/>Tru-Link Leadwire sets , 3-Lead , Grabber , IEC, 24in &amp; 2PIN</p>                |
| <p><b>ML5726</b><br/>One piece ECG cable with leadwires, 5-lead, IEC, snap, Round 17pin. Spacelab Medical, 90496 Ultraview, 90369 Ultraview, 90367 Ultraview ) Spacelabs: 90701, 90303, 90305, 90308, 90309, 90385</p>   |                        |    | <p><b>ML5727</b><br/>One piece ECG cable with leadwires, 5-lead, IEC, Grabber, Round 17pin. Spacelab</p> |
| <p><b>ML5728</b><br/>One piece ECG cable with leadwires, 3-lead, IEC, snap, Round 17pin. Spacelab Medical, 90496 Ultraview, 90369 Ultraview, 90367 Ultraview. Spacelabs: 90701, 90303, 90305, 90308, 90309, 90385</p>  |                      |   | <p><b>ML5729</b><br/>One piece ECG cable with leadwires, 3-lead, IEC, Grabber, Round 17pin. Spacelab</p> |
| <p><b>ML8708</b><br/>CABLES ECG / E.C.G.CABLES - E.C.G. ACCESSORIES - MULTIPARAMETER CABLES. Marque: PHILIPS - AGILENT - HEWLETT PACKARD (MONITORAGE) Type: VM 4 (VM4). 71782. CABLE MULTIPAR. PHILIPS-AGILENT. H.P.3 PARAM. ECG SPO2 Temp-IEC</p>                                   |   | <p>Cabo para Monitor Philips SureSigns VM4</p>  |  |
| <p><b>ML9055</b><br/>Cabo Dräger Précordial C/5 DERIVAÇÕES (P/monitor DRAGER REF. 5956367)</p>   |   | <p>Dräger 5-lead grabber set, InfinityR Telemetry, IEC1 European (European color code) Dräger Ref. 5956367</p>                    |  |
| <p>Schiller SCOPE CM8-TM7, 12 MALE PINS PROTECTED Schiller SCOPE CM8 – TM7<br/><b>See ML5221</b></p>   |   |    | <p>Schiller Scope CM8-TM7, 12 Male Pins Protected Schiller Scope CM8 – TM7<br/><b>See ML5542</b></p>     |
| <p><b>ML4408</b><br/>Cabo Monitor Mediana Hanox E-10302. One piece ECG cable with 3-Id leadwires, snap, AHA) (monobloco de 3 précordiais). Mediana Model YM-2002 ES Nellcor Hanox 010725 ref.SP193 *</p>   |   | <p>Cabo Monitor Mediana Hanox EM - 0209</p>  |                     |




|  |  |   |
|--|--|---|
| <p><b>ML2215</b><br/>Cables for mindray MEC 1000 (MEC1000) (3 LEADS MAIN CABLE WITHOUT LEADWIRE SPACELABS SCOPE 90701-90303-90305-90308-90309-90385 PROTECTED. ECG Patient Cables PRO 1000. 107327: 3-Lead ECG</p>                         |   | <p><b>ML2219</b><br/>Din safety ECG leadwires set,3-Lead ,Snap ,IEC ,40in ,DIN</p>   |
| <p><b>ML9209</b><br/>Cables for minday MEC 1000 (MEC1000) (5 LEADS MAIN CABLE WITHOUT LEADWIRE SPACELABS SCOPE 90701-90303-90305-90308-90309-90385 PROTECTED. ECG Patient Cables PRO 1000. 107327: 5-Lead ECG</p>                          |   | <p><b>ML2216</b><br/>Din safety ECG leadwires set,5-Lead ,Snap ,IEC ,40in ,DIN</p>   |
| <p><b>ML7094</b><br/>Lifelink safety ECG trunk cable, 5-lead, 9ft, 120° angle connector, 6PIN. Mindray, CSI, Datascope Spacelab, BCI, Nellcor, Welch Allyn, Critikon. See model Compatibilities Below ML5088 or ML8725</p>                 |   | <p><b>ML7095</b><br/>LL safety ECG leadwires set, 5-Lead , Clip , IEC ,40in , LL-2PIN</p>    |
| <p><b>ML7093</b><br/>Lifelink safty ECG trunk cable, 3-lead, 9ft, 120° angle connector, 6PIN. Compatible Brand Mindray, CSI, Datascope Spacelab, BCI, Nellcor, Welch Allyn, Critikon. See model Compatibilities Below ML5088 or ML8725</p> |  <p style="text-align: center;"><b>ML7093</b></p>  |  <p style="text-align: center;"><b>ML7092</b></p> <p><b>ML7091</b><br/>LL safety ECG leadwires set,3-Lead , Clip , AHA,40in ,LL-2PIN</p>  |
| <p><b>ML7089</b><br/>5-lead Trunk cable, AHA/ IEC, 6pin, without resistance. Compatible Brand Mindray, CSI, Datascope Spacelab, BCI, Nellcor, Welch Allyn, Critikon. See model Compatibilities Below ML8725 or ML5088</p>                  |   | <p><b>ML7079</b><br/>5-lead IEC leadwires, Snap</p>    |
| <p><b>ML8722</b><br/>Cables for minday MEC 1000 One piece ECG cable with 3-ld, TPU, IEC, Snap, 12ft &amp; 120° angle connector, 6pin. Cables for minday MEC 1000 Zoll 579-A-3ST. Patient Monitor CMS6000/CMS7000/CMS, EDAN</p>             |  <p style="text-align: center;"><b>ML8722</b></p>  <p style="text-align: center;"><b>ML4221</b></p> | <p><b>ML4221</b><br/>Cables for minday MEC 1000 One piece ECG cable with 3-ld, TPU, IEC, Clip, 12ft &amp; 180° straight connector .6pin Cables for minday MEC 1000 Zoll 579-A-3ST, EDAN</p>    |
| <p><b>ML8725</b><br/>Cables for minday MEC 1000 One piece ECG cable with 5-ld, TPU, IEC, Snap, 12ft &amp; 120° angle connector, 6pin. Cables for minday MEC 1000, Mek MP700. Zoll 579-A-3ST. EDAN Patient Monitor CMS6000/CMS7000/CMS</p>  |  <p style="text-align: center;"><b>ML8725</b></p>  <p style="text-align: center;"><b>ML5088</b></p> | <p><b>ML5088</b><br/>Cables for minday MEC 1000 One piece ECG cable with 5-ld, TPU, IEC, Clip, 12ft &amp; 120° angle connector. 6pin, Zoll 579-A-3ST. Patient Monitor CMS6000/CMS7000/CMS. EDAN</p>    |
| <p><b>ML5119</b><br/>Larsen &amp; Toubro (L&amp;T) 3-lead ECG cable with leadwire, 12pin , IEC , Snap</p>   |  <p style="text-align: center;"><b>ML5119</b></p>  <p style="text-align: center;"><b>ML8741</b></p> | <p><b>ML8741</b><br/>Corometrics 556, Critikon 8100T, Dinamap 8700, 9700 Series Lifepak 6, Nellcor N-200 N-250, N-1000, NPB-3900, NPB-4000, Nihon Kohden :BSM-8301A, BSM-8500A, BSM-8502A, OEC-6102A, 6105A, 7102A, 8108A, TEC-5200a, TEC-7100A/7200A/7300A</p>  |



|   |  |  |   |
|---|--|--|---|
| <p><b>ML1531</b><br/>Din safety ECG leadwires set,5-<br/>Lead , Clip , IEC ,40in ,DIN</p>   |     |   | <p><b>ML2213,</b><br/>Din safety ECG<br/>leadwires set,3-<br/>Lead ,Clip ,<br/>AHA ,40in ,DIN</p>   |
| <p><b>ML3304</b><br/>Cable for Fukuda Denshi FCP-2155.<br/>ECG Fukuda FX-2111, FX7202,<br/>ECG Fukuda Cardimax FX-7102<br/>Cardisunny Fukuda-ME X1000.<br/>ECG cables to FUKUDA DENSHI<br/>FUKUDA FCP-7101. Fukuda CP-104T<br/>DB15M&gt;<b>Banana 4.0</b></p> |     |   | <p><b>Banana 4.0</b><br/>FCP-2155.<br/>Cardimax<br/>FX-2111<br/>/FX-3010 /<br/>FX-4010,<br/>FCP-2155<br/>Cable CP-104T</p>                                |
| <p>SIRECUST 600 - SIRECUST<br/>732 ( p/ECG Siemens Sirecust<br/>1260) Sirecust 341-342 R<br/>Series 200-400-600-700-900-<br/>1200<br/>10 Male Contacts</p>  |     |   |   |
| <p>SIRECUST 600 - SIRECUST<br/>732 ( p/ECG Siemens Sirecust<br/>1260) Sirecust 341-342 R<br/>Series 200-400-600-700-900-<br/>1200<br/>10 Male Contacts</p>  |     | <p>Monobloc Cable</p>  |   |
| <p>AS73524MQ Stammkabel<br/>kompatibel zu Marquette Hellige . 5-<br/>adrig für Solar, Dash, Eagle, Länge<br/>2,75 m trunk cable compatible for<br/>Marquette Hellige, 5 leads, for Dash,<br/>Eagle, Solar, Length: 2,75 m</p>                                 |  | <p>AS414556-003<br/>Patientenleitung . Original<br/>Marquette Hellige, 5 adrig mit<br/>Klammer, Länge: 74 cm.<br/>Original Marquette Hellige<br/>patientleads 5 leads with<br/>clamp, length: 107 cm</p> | <p>Cabo ECG<br/>p/monitor<br/>EAGLE<br/>Mod. 3100</p>   |
| <p><b>ML1274</b><br/>ECG Cable MEK<br/>MP700, MP800<br/>One piece Cable with<br/>Leadwires, PAMO II<br/>3-Lead, Snap, <b>7pin AHA</b> ECG<br/>Cable M160246 M360697 AHA or<br/>IEC medical TPU cable</p>  |   |   | <p><b>ML1275</b><br/>ECG Cable<br/>MEK<br/>MP700, MP800<br/>One piece<br/>Cable with<br/>Leadwires,<br/>3-Lead, Snap,<br/><b>7pin IEC</b><br/>PAMO II</p> |
| <p><b>ML4999</b><br/>Philips Internal Defibrillation<br/>Paddles M4741A<br/>(989803109731). M4741A Extra<br/>Large 7.5 cm (2.95 in).<br/>7.5cm Switched Internal<br/>Paddles. Switched 7.5cm</p>  |   |   |   |
| <p><b>ML9056</b><br/>MS20094C7, Siemens SC5000,<br/>SC6000, SC7000, SC8000,<br/>SC9000 6-pod multimed<br/>Trunk cable, drager infinity</p>  |   |   |   |
| <p><b>ML9057</b><br/>MC093-LW-5IG, <b>5 Leads</b>,<br/>MS20094C7, 6-POD leadwire<br/>for MC093, 6-lead , Grabber ,<br/>IEC</p>  |   | <p><b>ML9057</b><br/>MC093-LW-5IG, <b>3 Leads</b>,<br/>MS20094C7, 6-POD<br/>leadwire for MC093, 6-lead<br/>, Grabber , IEC</p>   |   |



|   |  |  |   |
|---|--|--|---|
| <p><b>ML7705</b><br/>Cabo PC ref. 91-0056-00 p/ligação do PC ao aparelho de Monitorização da Tensão Arterial, Marca TRACKER NIBP REYNOLDS MEDICAL .</p>                     |     |  |   |
| <p><b>ML9805</b><br/>BPM-103 Biosys One-piece ECG cable, 3-Lead, Snap, AHA, 12ft, Redel-6PIN.<br/>Compatibilities; BIOSYS BPM-103 (OLD).</p>                                |     | <p><b>ML9806</b><br/>BPM-103 Biosys One-piece ECG cable, 5-Lead, Grabber, IEC, 12ft, Redel-6PIN.<br/>Compatibilities; BIOSYS BPM-103 (OLD).</p>                  |   |
| <p><b>ML1284</b><br/>Bionet BM3 Compact E, 8 pins. One piece cable with leadwires, 3-Lead monitoring IEC 8pin, Grabber.<br/>BIONET – BM3 Compact E. ECG-10033D</p>          |     |   |    |
| <p><b>ML8702</b><br/>Bionet one piece ECG cable with leadwire, 8pin Round Connector, 5-lead, IEC, Grabber.<br/>Bionet BM3 patient monitor, one piece Cable with 5-lead</p>  |     | <p><b>ML4989</b><br/>Bionet one piece ECG Cable with 3-lead AHA Snap lead wires. Bionet BM3 patient monitor, one piece Cable with 3-lead AHA Snap leadwires.</p> |    |
| <p><b>ML4808</b><br/>BIONET – BM3, Cabo ECG 3 Derivações P/Monitor Cardíaco BIONET – BM3. Bionet ECG trunk cable, 3-Lead, LL/3-&gt; Bionet 8pin. BIONET – BM3 Compact E</p> |    | <p><b>ML1809</b><br/>BIONET – BM3, Cabo ECG 3 Derivações P/Monitor Cardíaco BIONET – BM3. Bionet ECG trunk cable, 5-Lead, LL/3-&gt; Bionet 8pin.</p>             |   |
| <p><b>ML1459</b><br/>LL-22305-P, LL safety ECG leadwires set 5-Lead, Snap, AHA, 40in, LL-2PIN<br/>BIONET – BM3 Compact E</p>  |   | <p><b>ML4809</b><br/>LL-22363-P, LL safety ECG leadwires set 3-Lead, Snap, AHA, 40in, LL-2PIN. BIONET – BM3<br/>BIONET – BM3 Compact E</p>                       |   |
| <p><b>ML5606</b><br/>MEDTRONIC PHYSIO-CONTROL ECG cable, 5-Lead, Snap, IEC, 12ft, Round 12pin.</p>  |   |  |   |
| <p><b>ML5607</b><br/>M&amp;B One-piece Shield ECG cable, 3-Lead, Snap, AHA, 12ft, round 6PIN</p>  |   |   |   |
| <p><b>ML6701</b><br/>Hellige one piece ECG cable with leadwire, 5-lead, IEC, Grabber.</p>   |   | <p><b>ML6702</b><br/>Hellige one piece ECG cable with leadwire, 3lead, AHA, Grabber.</p>   |  |
| <p><b>ML6706</b><br/>Cardioserv Defibrillator One piece ECG cable with leadwire, 3-lead, IEC, Grabber. Cardioserv Defibrillator cardio serv</p>                             |  | <p><b>ML6706A</b><br/>Cardioserv Defibrillator One piece ECG cable with leadwire, 3-lead, IEC, SNAP. Cardioserv Defibrillator Cardio serv</p>                    | <p>Use with smu611/SM U612. Use with Hellige Cardioserv Defibrillator</p>             |



|  |   |  |   |
|--|---|--|---|
| <p><b>ML5520</b><br/>Nihon Kohden<br/>Lifescop OPV-1500<br/>Trunk cable, 3M 20pin -&gt;8pin socket. Nihon Kohden Lifescop OPV-1500</p>   |                          |   | <p><b>ML5519</b><br/>Nihon Kohden<br/>OPV-1500<br/>One piece ECG cable with Leadwire, 3-lead, IEC, Grabber, 3M 20pin.</p> |
| <p><b>ML1191</b><br/>Datex-Ohmeda AS/3, CS/3, Cardiocap , Saltite plus, S/5 light monitor, With Resistance Cabo ECG p/ monitor sinais vitais Datex Cardiocap II</p>  |                          | <p>one piece ECG cable with leadwire, 5-lead, IEC, Snap, 10 pin female 4, 7K resistance,</p>   |                                        |
| <p><b>ML1192</b><br/>Datex-Ohmeda AS/3, CS/3, Cardiocap , Saltite plus, S/5 light monitor, With Resistance Cabo ECG p/ monitor sinais vitais Datex Cardiocap II</p>  |                          | <p>one piece ECG cable with leadwire, 5-lead, IEC, Grabber, 10 pin female 4,7K resistance.</p>   |                                        |
| <p><b>ML1193</b><br/>Datex-Ohmeda AS/3, CS/3, Cardiocap , Saltite plus, S/5 light monitor, With Resistance Cabo ECG p/ monitor sinais vitais Datex Cardiocap II</p>  |                          | <p>one piece ECG cable with leadwire, 3-lead, IEC, Grabber, 10 pin female 4,7K resistance.</p>   |                                        |
| <p><b>ML1194</b><br/>Datex-Ohmeda AS/3, CS/3, Cardiocap , Saltite plus, S/5 light monitor, With Resistance Cabo ECG p/ monitor sinais vitais Datex Cardiocap II</p>  |                         | <p>one piece ECG cable with leadwire, 3-lead, IEC, SNAP, 10 pin female 4,7K resistance.</p>  |                                       |
| <p><b>ML5608</b><br/>Drager one-piece Shield ECG cable, 3-Lead, Snap, IEC, 12ft, round 10pin. Drager Medical, PM8010 , PM8014, Babyguard 8000 , Cato, Babylog 8000/8000 plus PB8800 Series 8000, Cicero, Parameterbox UM1/3/3.1 , RM1,</p> |                        |   |   |
| <p><b>ML6009</b><br/>Part n.13179, S&amp;W Artema One piece ECG cable, 3-Lead, Clip, 9ft &amp; 10PIN, IEC.<br/>Also see our code ML8511</p>  |                        | <p><b>ML6010</b><br/>S&amp;W (Artema) One piece ECG cable, 5-Lead, Clip, 9ft &amp; 10PIN, IEC.</p>  |                                      |
| <p><b>ML8533</b><br/>Part n.13179, S&amp;W Artema One piece ECG cable, 3-Lead, Snap, 9ft &amp; 10PIN, IEC.</p>   |                       | <p><b>ML6033</b><br/>S&amp;W (Artema) One piece ECG cable, 5-Lead, Snap, 9ft &amp; 10PIN, IEC.</p>   |                                      |
| <p><b>ML6755</b><br/>Medtronic 5433A<br/>Part nº 5433 A, Medtronic 5433A Pulse Generator Pacing Cable. Use with Medtronic 5318 and 5438 Pulse Generator Hyperstronic 2 pin to device.</p>  |                        |   | <p>Medtronic 5433 A. Pulse Generator Pacing Cable</p>   |
| <p>We have Discontinued Cables and spare parts</p>   | <p>Para qualquer tipo de Cabo de monitorização ECG, Sensores ou Componentes Agradecemos nos contactem</p> | <p>For other cables, Sensors or spare parts, please contact us</p>   |                                      |